

# MACROMAN™



*Unique Comfort  
& Controlled Pipetting*



## Stable, safe & easy to pick up

- > Integrated **bench stand** adjustable for **left and right** handed users.
- > **Keep your MACROMAN handy** and ready for use, without risk of contamination.

## Soft & fatigue free pipetting

- > **Soft-touch material** for maximum finger comfort.
- > **Large dosing button** to decrease pushing efforts.



# MACROMAN™

## Optimal comfort

**Integrated, adjustable stand** for easy handling when pipette attached

**Large & soft-touch** material dosing button

**Light** weight

Quiet operation & **vibration free**

**Ergonomic**-shape handle with grip

## Safety & easy maintenance

**All liquid contact parts are autoclavable**

**Easy-to-replace** filter

Replaceable batteries

## Performance

**For all plastic & glass pipettes** from 1 to 100mL

**Continuously adjustable speed:** 25mL in 4sec

Environmentally friendly NiMH battery  
for **8 hours operation**

**Precise & smooth dosing button control:** drop-by-drop, gravity mode  
and motorized speed.

## Convenient accessories

**Horizontal or vertical & magnetic charging stand** (with complete mounting kit) for busy laminar flow cabinet.

**Unique adaptor** for D5000 and D10mL DIAMOND® tips.



Charging stand



## MACROMAN

<b>Weight</b>	194g
<b>Battery</b>	2x 1.2V AAA NiMH, removable, 8-hour operation, rechargeable during use, discharged & charging light.
<b>Autoclavability</b>	Nose piece, Silicone pipette mount & filter rubber.
<b>Filter</b>	Replaceable, 0.45µm or 0.20µm standard & sterile.
<b>Pipette types</b>	Glass or plastic from 1mL to 100mL
<b>Power supply</b>	100-240V / 50-60Hz/ EU, US, UK, JAP or AUS plug version
<b>Certification</b>	CE, EMC, RoHS, WEEE
<b>Contents</b>	Complete MACROMAN + holder + 1 extra 0.45µm sterile filter + identification labels + EU, US, UK, JAP or AUS plug version

Description	Reference Number
MACROMAN, EU plug	F110751
MACROMAN, US plug	F110752
MACROMAN, UK plug	F110753
MACROMAN, JAP plug	F110754
MACROMAN, AUS plug	F110755
Charging stand	F1077501
0.45µm standard filter	F1077504
0.45µm sterile filter	F1077525
0.20µm standard filter	F1077503
0.20µm sterile filter	F1077526
Nose piece	F1077506
Filter rubber	F1077502
Pipette mount	F1077505
Adapter for D5000 and D10mL DIAMOND® Tips	F161287

## Distributed by

**Gilson provides 1-year warranty after the date of purchase**

[www.gilson.com](http://www.gilson.com) | [sales@gilson.com](mailto:sales@gilson.com) | [service@gilson.com](mailto:service@gilson.com) | [training@gilson.com](mailto:training@gilson.com)

**Gilson, Inc. World Headquarters** | 3000 Parmenter Street | P.O. Box 620027 | Middleton, WI 53562-0027 USA  
Tel: (1) 800-445-7661 or (1) 608-836-1551 | Fax: (1) 608-831-4451

**Gilson S.A.S.** | 19, avenue des Entrepreneurs | BP 145, F-95400 Villiers-le-Bel, France  
Tel: (33-1) 34 29 50 00 | Fax: (33-1) 34 29 50 20

**ISO9001 Certified**

LT800498E - © 2009 Gilson S.A.S. - All rights reserved - Printed in France, June 2009. Gilson, Inc. and Gilson S.A.S are ISO 9001 certified. Specifications subject to change without notice.

