

# HydraCare® MICRO-TOUCH®

## Micro-Touch® HydraCare®

Powder-free natural rubber latex  
non-sterile examination glove



## A unique combination of protection and skin care

To minimise the effects of continual glove use on hands, Ansell offers Micro-Touch® HydraCare®. This powder-free, latex examination glove contains an interior moisturiser that restores skin hydration and elasticity.

In addition, healthcare workers and patients alike appreciate the glove's pleasant scent. Along with low allergenicity, Micro-Touch® HydraCare® also offers high-quality protection, comfort and ease of use.

Every Healthcare Professional  
deserves **Optimal Protection**

**Ansell**

HydraCare®

# MICRO-TOUCH®

## Micro-Touch® HydraCare®

Powder-free natural rubber latex non-sterile examination glove

### A unique combination of protection and skin care

#### APPLICATIONS

All non-sterile medical applications requiring increased grip and skin moisturising for healthcare workers.

A pale green colour ensures easy product identification, while the glove's pleasant scent is much appreciated by healthcare workers and patients alike.

#### BUILT-IN SKIN CARE

To help maintain skin suppleness, the inner surface of the Micro-Touch® HydraCare® glove is treated with a glycerol moisturiser. This has a proven corrective effect on defective skin, restoring its hydration and elasticity.\*

#### SAFETY AND LOW ALLERGENICITY

Providing a secure barrier against viruses, micro-organisms, chemicals and detergents, Micro-Touch® HydraCare® meets or exceeds all international and national safety standards.

#### COMFORT, EASY DONNING AND SURE GRIP

With a length of 240 mm, the Micro-Touch® HydraCare® glove fits perfectly, due to its unique design. Single chlorination combined with silicone and the moisturiser guarantee easy donning, even with damp hands.

Extensive leaching phases and post-production washing cycles make sure the glove has a very low residual protein content. This, along with the fact that it is powder-free, significantly reduces the risk of Type I allergies. Furthermore, the low level of residual process chemicals and absence of thiurams minimise the risk of Type IV allergies.

## MICRO-TOUCH®

### Ansell Micro-Touch® examination gloves

The Micro-Touch® range of gloves includes versatile and innovative products covering all examination applications, with a solution for each specialty. Ansell has been especially attentive to customer requirements regarding barrier protection, comfort, guaranteed length, wide size range, raw material formulation and easy dispensing.

In the design and manufacture of the Micro-Touch® range, Ansell has drafted strict technical specifications. This ensures the high quality and consistency of all Micro-Touch® products. Proof of this can be found in the values for AQL, force at break, length and protein levels.

#### ANSELL MICRO-TOUCH®

Because every healthcare professional deserves optimal protection.

### Specifications

\* **Glycerol Accelerates Recovery of Barrier Function In Vivo** J. W. Fluhr, M. Gloor, L. Lehmann, S. Lazzarini, F. Distanto, E. Berardesca.

**Abstract:** Glycerol creates a stimulus for barrier repair and improves the stratum corneum hydration; stratum corneum hydration is not strictly related to barrier homeostasis and can be optimised by different mechanisms and pathways. The observed effects were based on the modulation of barrier repair and were not biased by the humectant effect of glycerol. As the glycerol-induced recovery of barrier function and stratum corneum hydration were observed even 7 days after the end of treatment, glycerol can be regarded as a barrier stabilising and moisturising compound.

Size	X-Small	Small	Medium	Large	X-Large
Code	553331	553332	553333	553334	553335

### Packaging

100 gloves per box,  
10 boxes per carton

X-Large : 80 gloves per box,  
10 boxes per carton



Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

Ansell Limited is a global leader in barrier protective products. With operations in the Americas, Europe and Asia, Ansell employs more than 11,000 people worldwide and holds leading positions in the natural latex and synthetic polymer glove and condom markets. Ansell operates in three main business segments: Occupational Healthcare, supplying hand protection to the industrial market; Professional Healthcare, supplying surgical and examination gloves to healthcare professionals; and Consumer Healthcare, supplying condoms and consumer hand protection. Information on Ansell and its products can be found at <http://www.ansell.eu>

**Ansell UK Ltd** 30-32 Mariner Lichfield Road Industrial Estate Tamworth B79 7UL  
United Kingdom Phone: +44 (0) 1564 711 802 Fax: +44 (0) 1564 711 344  
<http://www.ansell.eu> E-mail: [infouk@eu.ansell.com](mailto:infouk@eu.ansell.com)

**Ansell Healthcare Europe N.V.** Riverside Business Park Block J Boulevard International 55  
B-1070 Brussels Belgium Tel.: +32 (0) 2 528 74 00 Fax: +32 (0) 2 528 74 01  
Fax Customer Service: +32 (0) 2 528 74 03 <http://www.ansell.eu> E-mail: [info@ansell.eu](mailto:info@ansell.eu)

# Ansell