



## Prefilled Tubes and Bulk Beads: The Right Bead for Every Sample

Wide range of kits available
Variety of bead sizes and types to
meet your sample needs

TriplePure™ Quality
Nuclease, Protease and Nucleic
Acid free

Prefilled 2.0ml tubes or bulk beads
Convenience of prefilled tubes
Bulk beads for filling your own

Benchmark Scientific's homogenizer bead kits come in a variety of materials and sizes for processing a wide range of samples. Prefilled for convenience, the unskirted 2.0ml tube kits will fit in all popular bead homogenizers including Benchmark Scientific's BeadBug™ and BeadBalaster™ 24. Tested to ensure they are RNase, DNase, protease, and nucleic acid free, the TriplePure Zirconium Bead Kits are supplied in five bead sizes to homogenize everything from bacteria to large tissues. Acid washed silica and stainless steel bead kits as well as specialty bead mixes and bulk beads are also available.

<b>Bead material</b>	Bead Diameter	Prefilled 2.0ml tubes, pk/50	Bulk cat. no./bottle		
TriplePure Zirconium	0.1mm	D1032-01	D1132-01TP/250g		
TriplePure Zirconium	0.5mm	D1032-05	D11 <mark>32-05TP/250</mark> g		
TriplePure Zirconium	1.0mm	D1032-10	D1132-10TP/250g		
TriplePure Zirconium	1.5mm	D1032-15	D1132-15TP/250g		
TriplePure Zirconium	3.0mm	D1032-30	D1132-30TP/250g		
TriplePure Zirconium	0.1, 0.5, 1.0, 1.5, 3.0mm	D1032-SK*	NA		
Silica	0.1mm	D1031-01	D1131-01/200g		
Silica	0.5mm	D1031-05	D1131-05/200g		
Silica	1.0mm	D1031-10	D1131-10/200g		
Stainless Steel	2.8mm	D1033-28	D1133-28/1,000 balls		
Garnet, Zirconium	Garnet shards/6mm zirconium satellite	D1033-30G	NA		
Silica, Zirconium	Mixed size	D1034-MX	NA		
Empty tubes, with scre	D1031-T20				
Empty tubes, with scre	D1031-T21				
<b>BeadBug Bead Homog</b>	D1030				
BeadBlaster 24 Bead H	D2400				
*Starter kit - contains 10 or					

Benchmark Scientific PO Box 709 Edison, NJ 08818 Phone 908 769-5555 Fax 908 222-1864 email info@benchmarkscientific.com Visit our website: www.benchmarkscientific.com

	Bead material	Silica			Zirconium				Steel	Garnet with Zr Satellite	Silica and Zirconium	
	Bead diameter	0.1mm	0.5mm	1.0mm	0.1mm	0.5mm	1.0mm	1.5mm	3.0mm	2.8mm	0.5/6mm	Mixed
Animal	Adipose							Х				
	Artery							Х	Х			
	Bone									Х		
	Brain							Х	Х			
	Insect - fly								X	X		
	Insect - Tick								Х	Х		
	Insect - Mosquito						Х	Х	X			
	Gonad							Х	Х			
	Hair								X	Х	X	
	Heart							Х	х			
	Kidney							X				
	Liver											
	Lung							Х	Х			
	Muscle								Х	Х		
	Pancreas							Х	Х			
	Skin								Х		X	
	Spleen							Х	Х			
	Teeth								Х	Х		
Microbial	Algae		Х			Х						
	Bacteria	х	Х		X	Х						
	Fungi		Х	Х		Х	Х					
	Mold		Х			Х						
	Spores		Х	Х		Х	Х					
	Yeast		Х	Х		Х	Х					
Plant and Soil	Leaves								Х	Х		
	Nuts									Х		
	Plant, soft tissue								Х	Х		
	Plant, tough tissue								х	Х	Х	
	Seeds								Х	Х		
	Stems								х	Х		
	Soil							Х				Х
	Coral									Х		
E	Stool							Х				Х
Fecal	Scat								Х	Х		х

Above is a partial list of sample types and suggested beads. For additional information and technical support, contact Benchmark Scientific at info@BenchmarkScientific.com.



PO Box 709 Edison, NJ 08818 TEL: 908 769-5555 Fax: 908 222-1864 email: info@BenchmarkScietific.com Web: www.BenchmarkScientific.com