

To Whom It May Concern
Subject: SMARTLine™ (1062X) - EU Type Examination Module B

Date: 4th March 2021

Dear Madam, dear Sir,

Hereby, for SMARTLine™ Powder-free blue nitrile gloves, 240 (1062X series) we confirm we have initiated the following steps:

- Request to extend the existing EU Type Examination certificate to include our SMARTLine™ gloves, 1062X series.
- Preparation of the submission of the technical files and the product to be certified under the SHIELD Scientific umbrella.
- Upgrade the glove from Type C (K) to Type B (KPT) according to ISO 374-1:2016+A1:2018 based on the latest test results below.

| CAS N° | CHEMICAL | LETTER CODE | DR (%) | BTT (min) | PERFORMANCE LEVEL |
|-----------|---------------------------|-------------|--------|-----------|-------------------|
| 1310-73-2 | Sodium Hydroxide 40% | K | -5 | >480 min | Level 6 |
| 7722-84-1 | Hydrogen Peroxide 30% | P | 18 | >480 min | Level 6 |
| 50-00-0 | Formaldehyde Solution 37% | T | 34 | 34 min | Level 2 |

Note: Tested according to EN 16523-1:2015+A1:2018 and ISO 374-4:2019.

Should you have any further questions on the SHIELD Scientific disposable gloves, please do not hesitate to contact your SHIELD Scientific representative.

Trusting to have informed you sufficiently.

Best regards

SHIELD Scientific B.V.



Cisco Robles
General Manager

SHIELDskin™ and ORANGE NITRILE™ are trademarks of SHIELD Scientific.

©2007 Copyright SHIELD Scientific B.V. - All Rights reserved

SHIELD Scientific B.V. • Dr. Willem Dreeslaan 1 • 6721 ND Bennekom • The Netherlands - Phone +31 (0)317 700 202 •
Fax +31 (0)318 503 742 • E-mail: Info@shieldscientific.com • www.shieldscientific.com