



## ALL CAPILLARIES AND PISTONS (CPs) AVAILABLE IN TIPACK™!

Gilson is now offering CP25-50-250 in TIPACK packaging of 96 assembled and ready-to-use CPs.





## Save time and get even more protection using MICROMAN!

### READY-TO-USE

CP 25, 50 & 250 are now pre-assembled.

### EASY-TO-FIT


CPs come in 96 unit racks: fitting becomes straight forward!

### ERROR-FREE

Racks are color-coded for easy identification.

### ULTIMATE PROTECTION

No need to touch the CP anymore, just grab your MICROMAN and get started!



Full traceability thanks to the batch identification number



Error-free with color-coded racks

## microman®

The Best PIPETMAN® Partner For Problem Liquids and Hazardous Samples!



- ✓ Maximum accuracy and precision are guaranteed for any application involving pipetting viscous, volatile, dense or hazardous liquids such as radioactive samples.
- ✓ 100% protection against contaminants.

## Still, the genuine Gilson MICROMAN CPs!

### HIGH CHEMICAL RESISTANCE

Capillaries made of pure polypropylene.

### CLEVER DESIGN

Light and slim to fit even into the narrowest and deepest test tubes.

### TWO OFFERS TO SUIT YOUR NEEDS

CP 25, 50 & 250 available:

- In bulk version of 200 CP units to assemble.
- In TIPACK rack version of 96 assembled CPs (sold in packs of 10 TIPACKS).

## CP range

MICROMAN model	Related CP	Volume range (µL)	Ordering reference	Quantity per box	Available in TIPACK
M10	CP10	1 - 10	F148412	192	✓
	CP10		F148312	960	✓
	CP10ST*		F148413	192	✓
	CP10ST*		F148313	960	✓
M25	CP25	3 - 25	F148012	960	✓
			F148112	200	-
M50	CP50	20 - 50	F148013	960	✓
			F148113	200	-
M100	CP100	10 - 100	F148414	192	✓
	CP100		F148314	960	✓
	CP100ST*		F148415	192	✓
	CP100ST*		F148315	960	✓
M250	CP250	50 - 250	F148014	960	✓
			F148114	200	-
M1000	CP1000	100 - 1000	F148560	182	✓
	CP1000ST*		F148180	182	✓

\*ST = Sterilized

■ New reference



[www.gilson.com](http://www.gilson.com) | [sales@gilson.com](mailto:sales@gilson.com) | [service@gilson.com](mailto:service@gilson.com) | [training@gilson.com](mailto:training@gilson.com)

100% recycled paper

**Gilson, Inc. World Headquarters** | 3000 Parmenter Street | P.O. Box 620027  
Middleton, WI 53562-0027 USA  
Tel: (1) 800-445-7661 or (1) 608-836-1551 | Fax: (1) 608-831-4451

**Gilson S.A.S.** | 19, avenue des Entrepreneurs | BP 145, F-95400 Villiers-le-Bel, FRANCE  
Tel: (33-1) 34 29 50 00 | Fax: (33-1) 34 29 50 20



**ISO9001 Certified**

LT800502 © 12/2009 Gilson SAS.  
All rights reserved. All trademarks are the property of Gilson Inc. and its subsidiaries.