



ALL CAPILLARIES AND PISTONS (CPs) AVAILABLE IN TIPACK™!

Gilson is now offering CP25-50-250 in TIPACK packaging of 96 assembled and ready-to-use CPs.





Save time and get even more protection using MICROMAN!

READY-TO-USE

CP 25, 50 & 250 are now pre-assembled.

EASY-TO-FIT


CPs come in 96 unit racks: fitting becomes straight forward!

ERROR-FREE

Racks are color-coded for easy identification.

ULTIMATE PROTECTION

No need to touch the CP anymore, just grab your MICROMAN and get started!



Full traceability thanks to the batch identification number



Error-free with color-coded racks

microman®

The Best PIPETMAN® Partner For Problem Liquids and Hazardous Samples!



- ✓ Maximum accuracy and precision are guaranteed for any application involving pipetting viscous, volatile, dense or hazardous liquids such as radioactive samples.
- ✓ 100% protection against contaminants.

Still, the genuine Gilson MICROMAN CPs!

HIGH CHEMICAL RESISTANCE

Capillaries made of pure polypropylene.

CLEVER DESIGN

Light and slim to fit even into the narrowest and deepest test tubes.

TWO OFFERS TO SUIT YOUR NEEDS

CP 25, 50 & 250 available:

- In bulk version of 200 CP units to assemble.
- In TIPACK rack version of 96 assembled CPs (sold in packs of 10 TIPACKS).

CP range

| MICROMAN model | Related CP | Volume range (µL) | Ordering reference | Quantity per box | Available in TIPACK |
|----------------|------------|-------------------|--------------------|------------------|---------------------|
| M10 | CP10 | 1 - 10 | F148412 | 192 | ✓ |
| | CP10 | | F148312 | 960 | ✓ |
| | CP10ST* | | F148413 | 192 | ✓ |
| | CP10ST* | | F148313 | 960 | ✓ |
| M25 | CP25 | 3 - 25 | F148012 | 960 | ✓ |
| | | | F148112 | 200 | - |
| M50 | CP50 | 20 - 50 | F148013 | 960 | ✓ |
| | | | F148113 | 200 | - |
| M100 | CP100 | 10 - 100 | F148414 | 192 | ✓ |
| | CP100 | | F148314 | 960 | ✓ |
| | CP100ST* | | F148415 | 192 | ✓ |
| | CP100ST* | | F148315 | 960 | ✓ |
| M250 | CP250 | 50 - 250 | F148014 | 960 | ✓ |
| | | | F148114 | 200 | - |
| M1000 | CP1000 | 100 - 1000 | F148560 | 182 | ✓ |
| | CP1000ST* | | F148180 | 182 | ✓ |

*ST = Sterilized

■ New reference



www.gilson.com | sales@gilson.com | service@gilson.com | training@gilson.com

100% recycled paper

Gilson, Inc. World Headquarters | 3000 Parmenter Street | P.O. Box 620027
Middleton, WI 53562-0027 USA
Tel: (1) 800-445-7661 or (1) 608-836-1551 | Fax: (1) 608-831-4451

Gilson S.A.S. | 19, avenue des Entrepreneurs | BP 145, F-95400 Villiers-le-Bel, FRANCE
Tel: (33-1) 34 29 50 00 | Fax: (33-1) 34 29 50 20



ISO9001 Certified

LT800502 © 12/2009 Gilson SAS.
All rights reserved. All trademarks are the property of Gilson Inc. and its subsidiaries.