



FAST MANUAL MICROPLATE PIPETTING



Save Time And Increase Your Productivity

With minimal setup you can fill 96- or 384- well plates in less than one minute by reducing the number of manipulations for filling a microplate.

Simplifies Sample Prep

PLATEMASTER® requires no programming, no power and minimal training, so you can start pipetting right away.

Convenient For Any Microplate Application

From plate transfers to supernatant, achieve accuracy and precision on the bench or under the hood.

PLATEMASTER is an amazingly easy-to-use, accurate solution for high throughput pipetting of 96- and 384-well microplates.

| PLATEMASTER | |
|----------------------------|---|
| PLATEMASTER P20 (F110761) | 0.5 µL - 20 µL Adjustable in 0.1 µL increments |
| PLATEMASTER P220 (F110762) | 2 µL - 220 µL Adjustable in 1 µL increments |
| Dimensions | 300 mm x 450 mm x 400 mm (W x L x H) |
| Weight | 16 kg |

| PLATEMASTER ACCESSORIES | |
|-------------------------|--|
| Description | Part Number |
| F1077602 | 384-well adapter, positioned done by wheel. |
| F1077603 | 384-well adapter, positioned done by hand. |
| F1077604 | Alu-heater block for PLATEMASTER®, 96 x 0.2 mL, for PCR tubes, stripes & PCR plates Ø 7.8 mm |
| F1077605 | Pipetting head height adjusters (2 units) |
| F1077606 | Lubrication box of O-rings PLATEMASTER® |

CARROUSEL™ Pipette Stand

- Protects pipettes from shocks and contamination.
- Holds up to seven Gilson pipettes for space saving.



Power CARROUSEL™

- PIPETMAN M charging stand.
- Charge up to 5 at once in one hour.
- Keeps your PIPETMAN M always operational.



Single™ Pipette Holder

- Holds your single or multichannel PIPETMAN.
- Sticks to any smooth, dry surface.
- Prevents liquid from running inside shaft.



Trio™ Pipette Stand

- Safely stores and protects three PIPETMAN.
- Helps organise your workbench.



Universal Multichannel Stand

- Holds any brand of 8 or 12 channel pipette securely.
- Compact design for space savings.



Gilson Reagent Reservoirs

- Disposable, ideal for dispensing reagent with multichannel pipettes.
- 25 mL and 50 mL available.



PIPETTE ACCESSORIES

PIPETMAN COMFORT HANDLE

- Clips on to PIPETMAN Classic, Neo and G models.
- Reduces the force needed to grip the pipette for maximum comfort.



Coloris™ Identification Clips

- Identify or personalise your PIPETMAN.
- Write your name or application to avoid mix-ups and then snap into place.



JIMMY™ Microtube Opener

- Clips on to pipette.
- Increases pipetting efficiency.
- Lessens the fatigue associated with manual tube opening.



PIPETMAN ACCESSORIES

| Description | Part Number |
|----------------------------------|-------------|
| Carrousel, stand for 7 pipettes | F161401 |
| Power Carrousel | FB1001 |
| Single pipette holder | F161406 |
| Trio, stand for 3 pipettes | F161405 |
| Universal Multichannel Stand | F161417 |
| Reagent Reservoir 25 mL, 100 | F267660 |
| Reagent Reservoir 50 mL, 100 | F267670 |
| COMFORT HANDLE pack of 3, blue | F161486 |
| COMFORT HANDLE pack of 50, blue | F161485 |
| COMFORT HANDLE pack of 3, white | F161488 |
| COMFORT HANDLE pack of 50, white | F161487 |
| COMFORT HANDLE stickers, 120 | F161499 |
| 10 Coloris clips, mixed colours | F161301 |
| 10 Coloris clips, red | F161302 |
| 10 Coloris clips, yellow | F161303 |
| 10 Coloris clips, green | F161304 |
| 10 Coloris clips, blue | F161305 |
| 10 Coloris clips, white | F161306 |
| The JIMMY Microtube Opener | F144983 |