



INCREASED COMFORT AND PRODUCTIVITY

Air-Displacement Pipettes





Adjustable Tip Ejector

Designed for performance and comfort, PIPETMAN® L offers a comprehensive range of models, adapted to all your needs.

Convenient Use

Lightweight and balanced, PIPETMAN L combines an ergonomic design with minimized pipetting forces, giving you the ability to pipette longer without fatigue.

Right-handed or left-handed? No problem. Adjust the location of the tip ejector to the most comfortable position for you.

Secure Pipetting

Reduce pipetting errors and avoid accidental volume changes during pipetting cycles with a true volume-locking mechanism for improved accuracy and precision.

PIPETMAN L variable volume models are all fully autoclavable without disassembly for maximum convenience and protection from cross-contamination.

Extended Range

Available as a fixed or variable volume, the PIPETMAN L line covers a full volume range from 0.2 μ L to 10 mL for single-channel models, and from 0.5 μ L to 1200 μ L for multichannel models. Precise, accurate, and reliable, they can be intensively used with confidence in many applications.

Four PIPETMAN L multichannel models (product description ending with VR), covering a volume range from 20 μ L to 300 μ L, are equipped with V-rings at the bottom of the tip holders for a proper fit with most tip brands.





Choose your Tip Ejector

For single channel models, choose between a model with a stainlesssteel tip ejector or a plastic ejector depending on your needs and applications.*

Large volume models are equipped with a plastic tip ejector to avoid any direct contact with the tip.



Choose your Model

Choose your channel number and volume range, up to 1200 μ L for various protocols.

Four PIPETMAN L multichannel models equipped with V-rings (VR) offer an improved compatibility with tips from other brands.





Choose your Pipette Stand

Different pipette stands and holders allow you to store pipettes vertically, protecting them from contamination.









PIPETMAN L single channel models are available in kits, offering the essential pipetting tools to start an experiment right away.



For routine testing applications, PIPETMAN L is also available as a fixed volume version, covering a volume range from 1 μ L to 5 mL.



PIPETMAN DIAMOND® Tips are designed to fit your PIPETMAN perfectly and guarantee maximum accuracy and precision.

For ordering information, visit www.gilson.com

PIPETMAN L Single Channel — Variable Volume Models						
Model	PIPETMAN® DIAMOND Tips		Part Number	Volume (μL)		
P2L	D10 DL10	DF10ST DFL10ST	FA10001P	0.2 - 2		
			FA10001M			
P10L	D10 DF10ST DFL10T	DF10ST	FA10002P	0.5 - 10		
PIOL		FA10002M	0.5 - 10			
P20L	D200	DF30ST	FA10003P	2 - 20		
PZUL	D200	DF3051	FA10003M			
P100L	D200 DF100ST	DE100CT	FA10004P	10 - 100		
PIOOL		DFI0031	FA10004M			
P200L	D200 DF200ST	FA10005P	20 - 200			
PZOOL	D300	DF300ST	FA10005M	20 - 200		
P1000L	D1000 D1200	DF1000ST DF1200ST	FA10006P	100 - 1000		
			FA10006M			
P5000L	D5000		FA10007	500 - 5000		
P10mLL	D10mL		FA10008	1000 - 10000		

For a pipette with a plastic tip ejector, the part number ends with the letter P. For a pipette with a stainless steel tip ejector, the part number ends with the letter M, from P2L to P1000L.

PIPETMAN L Multichannel Models					
Model	PIPETMAN® DIAMOND Tips		Part Number	Volume (μL)	
P8x10L	D10 DL10	DF10ST DFL10ST	FA10013	0.5 - 10	
P12x10L			FA10014		
P8x20L	DL10 D200	DFL10ST DF30ST	FA10009	2 - 20	
P12x20L			FA10010		
P8x200L	D200 D300	DF200ST DF300ST	FA10011	20 - 200	
P12x200L			FA10012		
P8x200LVR	D200 D300	DF200ST DF300ST	FA10035	20 - 200	
P12x200LVR			FA10036		
P8x300L	D200	DF200ST	FA10015	20 - 300	
P12x300L	D300	DF300ST	FA10016		
P8x300LVR	D200 D300	DF200ST DF300ST	FA10037	20 - 300	
P12x300LVR			FA10038		
P8x1200L	D1200	DF1200ST	FA10039	100 - 1200	
P12x1200L	D1200		FA10040		

Pipette Stands				
Description	Part Number			
CARROUSEL Pipette Stand (7 pipettes)	F161401			
TRIO Pipette Stand (3 pipettes)	F161405			
Universal Multichannel Stand	F161417			
SINGLE Pipette Holder	F161406			

