

SAPPHIRE Pipette Tips (Bulkware)

Demands are constantly changing. For this reason, we have completely revised and further developed our pipette tip range. We are now pleased to present our new, high-quality and clearly structured pipette tip portfolio.

With their innovative design and perfect fit, Sapphire pipette tips stand for the best possible user comfort and optimal results in everyday laboratory work.

Due to the high demand, we now offer our new tips in advance as bulkware packed in user-friendly bags. Corresponding empty racks are available separately.

The entire portfolio with many other variants, for example with filter, as low retention, in racks or as refill units, will be launched in the course of the next few months.





KEY FACTS

- Universal fit for all common pipettes
- / Easy-to-read graduation for quick volume checks
- / Made from highly transparent medical grade polypropylene
- Slim tip shape reduces contamination risk when pipetting
- / Thin-walled design for a perfect fit and reduced loading and ejection forces



SAPPHIRE PIPETTE TIPS (BULKWARE)



The bulk packaging of our new Sapphire pipette tips stands for maximum user comfort in daily laboratory use.

Arguments that convince:

- / High-quality and easy-to-open bags for convenient, hygienic removal of the pipette tips
- / Reclosability: after opening, the bags can be securely closed again
- / Colour-coded bag labelling for quick volume identification
- / Fully autoclavable

Sapphire Pipette Tips (Bulkware)

Tip colour: natural, Graduation: yes, non-sterile

Item No.	Nominal volume	Volume range	Matching rack	Rack insert colour	Packaging option	Qty. inner/outer
771350	10 µl	0.2 - 10 μΙ	770310	● red	bag	1,000 / 10,000
772350	10 µl extended	0.5 - 10 μΙ	770310	• red	bag	1,000 / 5,000
775350	200 μΙ	5 - 200 μΙ	770320	yellow	bag	1,000 / 15,000
776350	300 μΙ	10 - 300 μΙ	770330	green	bag	1,000 / 10,000
777350	1000 μΙ	50 - 1000 μΙ	770340	blue	bag	1,000 / 5,000

SAPPHIRE RACKS (UNFILLED)



To match our new pipette tips in bags, we offer transparent, unfilled racks with coloured rack inserts for easy identification of the tip volume. Your advantages:

- / Autoclavable and reusable
- / Easy handling due to optimal opening and closing
- / For use with single-channel and multi-channel pipettes
- / Stackable and automation-friendly

Sapphire Racks for Pipette Tips (unfilled)

Rack colour: natural, with lid and coloured rack insert, non-sterile, for 96 tips

Item No.	Suitable for	Rack insert colour	Qty. inner/outer
770310	10 μl / 10 μl extended	● red	10 / 60
770320	200 μΙ	yellow	10 / 60
770330	300 µl	green	10 / 60
770340	1000 μΙ	blue	10 / 40

HOW TO CHOOSE YOUR NEW SAPPHIRE BULKWARE TIPS

		PREVIOUS	NEW		
		Greiner Bio-One & Sapphire Tips	Sapphire Tips		
Volume	Packaging option	Item No.	Item No.		
10 µL	bag, non-sterile	771290 771291	771350		
10 µL EXTENDED	bag, non-sterile	765290 771250	772350		
200 µL	bag, non-sterile	685290 685295 739290 739291 739296	775350		
250 μL	bag, non-sterile	739250	→ see 300 µl		
300 μL	bag, non-sterile	738250	776350		
1000 μL	bag, non-sterile	686290 686295 740290 740291 740296	777350		



