

# OMNI Tips™ plastic disposable generator probes

The OMNI Tips' patented design eliminates the hassle of disassembly and maintenance of traditional stainless steel generators.

Proprietary OMNI Tips™ Probes are available in both soft tissue (deaggregation) and hard tissue (frozen) versions. Engineered from durable plastics, OMNI Tips™ are reusable (autoclavable) and are so economical they can be disposed of after each use to eliminate cross contamination.

The outer tube of OMNI Tips™ is made of polycarbonate. The inner shafts are made of Ultem 1000, a polyetherimide. The OMNI Tips' patented design eliminates the hassle of disassembly and maintenance of traditional stainless steel generators. Treatable volume - 250 µL to 30 mL.

## Features

- **Economical:** <1% of the cost of stainless steel probes
- **Eliminates Cross Contamination**
- **Ideal for PCR, RNA, and DNA sample preparation**
- **Maintenance Free**
- **Sterile OMNI Tips™ are Electron Beam treated to ensure they are free of contaminants.**

*"With OMNI Tips™ I can be certain that all the tissue sample has been homogenized; there is no sample loss. I found that the OMNI Tips™ hard tissue version worked as well as the stainless steel probe."*

— Rob T., National Institutes of Health

## NEW Sterile OMNI Tips™ now available

### A selection of our best-selling plastic probes:

#### SOFT TISSUE PLASTIC GENERATOR PROBES

	ORDER NO.
OMNI Tips™ Plastic Generator Probes, Soft Tissue 25 pack	30750
OMNI Tips™ Plastic Generator Probes, Soft Tissue 50 pack	32750
OMNI Tips™ Plastic Generator Probes, Soft Tissue 100 pack	34750

#### HARD TISSUE PLASTIC GENERATOR PROBES

	ORDER NO.
OMNI Tips™ Plastic Generator Probes, Hard Tissue 25 pack	30750H
OMNI Tips™ Plastic Generator Probes, Hard Tissue 50 pack	32750H
OMNI Tips™ Plastic Generator Probes, Hard Tissue 100 pack	34750H
<b>Sterile OMNI Tips™ Generator Probes, Hard Tissue 25 pack</b>	<b>30750H-S</b>

**Note:** A1000SB adapter is required for OMNI Tip™ use on OMNI GLH, GLH 850, PDH and TH (with black nose ring).

**Note:** 12-120 OMNI Tip™ adapter is required for OMNI Tip™ use on the OMNI TH115 and TH220 (with blue nose ring).



**OMNI**  
International

935 Cobb Place Blvd  
Kennesaw, GA 30144  
Phone: 770.421.0058 • Fax: 770.421.0206  
[www.omni-inc.com](http://www.omni-inc.com)

Rev. F