

Product name : <b>Vacucap™</b>	Catalogue no.: <b>T402-16 series</b>
Edition date : 2004-02-06	Revision : 2013-09-25
Approved by : Annette Roy, QA coordinator	
<i>This document replaces any previous version</i>	

## 🔗 **Product identification**

- Designation: **Vacucap™ tubes closures**

## 🔗 **Production Information**

- Producer's name: **Simport Scientific**  
➤ Address: 2588 Bernard Pilon, Beloeil (Québec), J3G 4S5  
➤ Tel. no / web site: (450) 464-1723 / [www.Simport.com](http://www.Simport.com)

## 🔗 **Packaging**

- Case Cardboard box containing 6 bags of 1000 units, for a total of 6000 caps

## 🔗 **Traceability.**

- Lot no  
➤ Composition: 8 digits  
➤ Location: ① On exterior case label.  
(The lot number can be found in different locations) ② On inserted label inside case.  
③ On inner package.

## 🔗 **Standard conformity.**

- **MIL-STD-105E** Sampling and inspection procedures  
➤ **FDA** Resin conforms to FDA 21 CFR 177.1520  
➤ **Coneg** Plastics and colorants are in conformity with CONEG standards

## Quality assurance

---

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements
- Closure verification

## Specifications

---

- Disposable caps for 16mm glass test tubes and plastic culture tubes.
- Material: Low Density polyethylene
- Temperature range: -80 °C to 100 °C
- Non autoclavable

## Storage conditions

---

- Store in cool dry place.
- Avoid temperature variations and humidity.
- Protect from any possible contamination.

## Recommended use

---

- Check cap closure when using biohazard material and/or chemical reagents.
- Follow chemical resistance chart recommendations
- Follow manufacturers instruction for automated equipment

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*