

► Thermo Scientific Sterilin Polystyrene Containers, 7mL Bijou



CE

Aseptic

Centrifugation

Sterilin Polystyrene Containers are ideal for small volume samples.

details

- Leak-tested in accordance with BS EN 14254 Annexe D
- Available pre-filled with boric acid (0.09g) for the preservation of urine samples
- Available with removable glass coverslip — ideal for work on the culture of intracellular organisms (e.g., Chlamydia trachomatis, Rickettsia, etc.)
- Suitable for centrifugation at 7,200 x g
- Aseptically manufactured
- *In vitro* use only

Ordering Information: Not available in all geographic regions

Certifications: CE marked in accordance with the European Directive 98/79/EC

Sterilin Polystyrene Containers, 7mL Bijou

Cat. No.	Description	Label	Material Base/Cap	O.D. x H	Sterility	Case Qty.
129A	Bijou	None	PS/PE	22.5 x 50.2mm	AS	700
129B	Bijou	Printed	PS/PE	22.5 x 50.2mm	AS	700
129BBAC	7mL Bijou + boric acid	Printed	PS/PE	22.5 x 50.2mm	AS	700
129AX/1	Bijou + coverslip	None	PS/PE	22.5 x 50.2mm	AS	700