

## Technical Data Sheet

**Product:** Steril stopper

**Code-no.:** 890 55 42

---

**Material:** Silicosen®

**Quantity of airflow:** 30 – 120 ml/ min

**Evaporation of Moisture:** 0,25 ml in 10 days at 30°C

**Sterilisation:** Dry sterilized (at 180°C)  
and high-pressure steam sterilized

**Cleaning:** Methods including rough washing in cold or hot water followed by boiling or ultrasound washing. Also neutral detergent is the suitable detergent. But note that all detergent must be completely rinsed off.

Hirschmann Laborgeräte  
GmbH & Co. KG



.....  
Roland Benz  
Quality Manager

Eberstadt, 21.02.2020