

## Technical Data Sheet

**Product:** Steril stopper Bio- Silico  
**Code-no.:** 890 59 ..

---

**Material:** Bio-Silico

**Quantity of airflow:** 3.600 – 4.500 ml/ min

**Evaporation of Moisture:** 0,5 ml in 10 days at 30°C

**Sterilization:** High-pressure steam sterilization  
(121 °C, 20 min.)  
Dry air sterilization (120°C, 20 – 30 min.)

Hirschmann Laborgeräte  
GmbH & Co. KG



.....  
Roland Benz  
Quality Manager

Eberstadt, 03.09.2020