



**Date: 01/01/2024**

**Re: Expiration Dates / Sterilization - Sterileware® Product Line**

Bel-Art products do not carry expiration dates.

The products remain sterile as long as the barrier remains unbroken preventing bacteria and or other airborne particulates from entering into the sterile item.

It is recommended that the items be properly stored to maintain the packages integrity.

Items should be stored in a clean and dry (non-humid) environment at ambient temperature while protected from direct sunlight or harsh chemicals until ready for use.

It is also recommended that the principal processing technique of first in – first out be used while maintaining inventory. Any definitive shelf life should be governed by the end user’s protocol.

Devices are generally regarded as sterile as long as package integrity has not been compromised. Also notable is that the manufacturing and sterilization processes have remained the same since the line’s introduction over 10 years ago.

I certify that I am authorized to declare that the statement made with this Certificate of Compliance is true and accurate.

I certify that I am authorized to declare that this statement is true and accurate.

A handwritten signature in blue ink, appearing to be 'Nancy Trogani', written over a horizontal line.

Signature: \_\_\_\_\_

***Nancy Trogani***

Technical Support Analyst

[Bel-Art](#) | [Wilmad](#) | [Ableware](#)

502 Valley Road Suite 201, Wayne, NJ 07470 USA

Tel: +1-973.694.0500 | U.S. +1-800.443.4926 | [info@belart.com](mailto:info@belart.com)

**Corporate Office:** 935 Mearns Road, Warminster, PA 18974 USA

[scientificproducts.com](http://scientificproducts.com)

**Bel-Art | FTS | Genevac | Hotpack | Hull i-Dositecno | Service | Training | VirTis | Wilmad**

