


CATALOG NO.	*3005-42/-70	
SIZE (NOM VOL. (mL))....	500	
VOL. AT NECK CAP. (mL)...	670	
MATERIALS	BOTTLE.....	PETG
	CLOSURE.....	WHITE PP
	CLOSURE LINER..	LDPE FOAM
NECK FINISH.....	38-430	
THREAD (T.P.I.).....	6	
NOM. WEIGHT W/ CLOSURE (gms).....	85	
	INCH	CM
A (NECK I.D.).....	1.11	2.8
B (HT. W/O CLOSURE).....	6.80	17.2
C (HT. W/ CLOSURE).....	6.94	17.6
D (SQUARE).....	3.0	7.7
W.T. (NOM. WALL THK).....	.020	.05

* 3005-42 IS LAB-PACK, STERILE
3005-70 IS BULK-PACK, STERILE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC-2	APP'D	PVC	PSS 3005-XX	REV.
CK'D	PVC		1/9/06		0106
DESCRIPTION: BIOTAINER BOTTLE, PETG, 500ml					