


CATALOG NO.	* 3230-20/-42	
SIZE NOM VOL. (mL)	2000	
NOMINAL BRIM CAP. (mL)	2,475	
MATERIALS	BOTTLE	PETG
	CLOSURE	WHITE PP
	CLOSURE LINER	LOPE FOAM
NECK FINISH	48mm DOUBLE START TH'D	
THREAD (T.P.I.)	3	
NOM. WEIGHT W/ CLOSURE (gms)	220	
	INCH	CM
A (NECK I.D. AT TOP)	1.44	3.6
B (HT. W/O CLOSURE)	10.31	26.2
C (HT. W/ CLOSURE)	10.45	26.5
D (SQUARE)	4.55	11.6
W.T. (NOM. WALL THK)04	0.1

* 3230-20 IS BULK-PACK, STERILE
3230-42 IS LAB-PACK, STERILE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC ²	APP'D	PVC	PSS 3230-XX	REV. 0106
CHK'D	PVC		1/9/06		
DESCRIPTION:					
BIOTAINER BOTTLE, PETG, 2L					