

<Specifications>

• • Electrophoresis Cell

Overall dimensions	183mm(W)x59mm(H)x162mm(L)
Material	PPHOX
Solution volume	Approx. 270-320 ml (includes gel tray and gel)
Multi-sample	Multi-channel pipette compatible

• • Safety Lid

Overall dimensions	197mm(W)x38mm(H)x170mm(L)
Safety interlock system	Without the lid, main power can not be operated.

• • Power Supply

Overall dimensions	75mm(W) x 62mm(H) x 170mm(L)
Weight	410g
Input voltage	AC100-240V (Internationally compatible), 50/60 Hz
Output voltage	135V 100V 50V 25V 70V 35V 18V Constant peak voltage of 140V and duty control
Timer	Can be set between 0 and 99 min. and continuous mode. Temporary shutdown supported
Memory function	Automatic memory (the last used V & T)

• • One-Gel Trays

Overall dimensions of small gel tray	130mm(W) x 16.5mm(H) x 59.5mm(L)
Overall dimensions of large gel tray	130mm(W) x 24mm(H) x 122mm(L)
Quantity	2 small gel trays; 1 large gel tray
Material	Heat Resistance Materials(For 100• Agarose)

• • One-Comb-HR

Multiple number of wells	13 - 26 wells (13wells: 9mm spacing, 26wells: 4.5mm spacing)
Quantity	4 (supports making 13 or 26 wells)
Material	Heat Resistance Materials(For 100• Agarose)

• • One-Gel Casting Stand-HR Makes 2 small gel trays or 1 large gel tray

Overall dimensions	140mm(W) x 20mm(H) x 125mm(L)
Material	Heat Resistance Materials(For 100• Agarose)

• • This product is subject to change without notice. • • "Mupid" is a registered trademark of ADVANCE Co., LTD.



IMPORTANT

To ensure safety, make sure you have thoroughly read the Instruction Manual before use.