



4 Reihen Noppen je Teilkreis
4 rows of knobs per pitch circle

Gewinde/Thread:
DIN 168-GL 45x4, 3 Gang/turn
0303803-GZA-01694--..

- 900
- 800
- 700
- 600
- 500
- 400
- 300
- 200
- 100



DURAN®

1000 ml



Made in Germany
DWK-LifeSciences.com

Retrace Code
00000000

Color: white
Farbe: weiß
0371401-SZ1-A

(22)

Dimensions not fitted with specific tolerances are to be seen as geometric dimensions only and are not subject to any tests/checks during production. These dimensions are approximate values only and could vary in the final products.

The reproduction, distribution and utilization of this document as well as the communication of its contents to others without express authorization is prohibited. Offenders will be held liable for the payment of damages. All rights reserved in the event of the grant of a patent, utility model or design.

21 801 54 14	Komplett PBT Kappe und ETFE Ring / with PBT Cap and ETFE Ring		
112 57 07	mit Staubschutzkappe / with Dust Cap		
21 801 54 06	Laborflasche 1000ml GL45 / Laboratory Bottle 1000ml GL45		
21 801 54 55	Komplett PP Kappe und Ring / with PP Cap and Ring		
21 806 54 02	braun / Amber		
21 806 54 51	braun; komplett / Amber; with PP Cap and Ring		
Art.-Nr./Art.-No.	Benennung/Title		
P	Logotausch DWK	14-MAR-18	Becker
N	Art.-Nr. 21 801 54 14 ergänzt; diverse Art.-Nr. entfernt	23-DEC-2016	Geier
Index	Aenderung / changes	Datum / date	Name
zul. Abweichung allowed dev.	Werkstoff / material	Stück quantity	Gewicht weight
Oberfläche surface	DURAN		Pos. position
	Benennung / title	Maßstab scale	
bearb. drawn	Datum / date	Name	Laborflasche mit GL45 Gewinde; 1000ml
geprüft checked	23.12.2016	Geier	Laboratory Bottle 1000ml, GL45
	Art.-Nr. see Tab.		1 : 1
	Form-Nr. mould number	1694/.1L	A3
	Zeichnungsnummer / drawing number	0303903-GZA-01694	Revision Blatt/sheet
		P	1/1
	Ers. / repl.f.	Ers. / repl.d.	
	Schutzvermerk ISO 16016 beachten / Refer to protection notice ISO 16016		INV

DWK
LIFE SCIENCES
DURAN Produktions GmbH & Co. KG
www.DWK.com