

THE LIGHT PROTECTED: DURAN® LABORATORY GLASS BOTTLE IN AMBER

07

- The amber DURAN® laboratory glass bottle also protects media from light radiation with a wavelength between 300 and 500 nm. It is, therefore, particularly suitable for use with light-sensitive media and for the long-term storage of substances.
- UV protection up to 500 nm wavelength
 - Consistent external application of the colour: internal surfaces have standard DURAN® glass properties
 - Uniform amber colour as a result of innovative technology
 - Special thermal method also results in long durability and good chemical resistance of the amber colour
 - Retrace code
 - Autodisable
 - Sizes: 10 ml to 20,000 ml

Order No.	Capacity (ml)	DIN-Thread (GL)	d (mm)	h (mm)	Pack Quantity
21 806 08	10	25	36	50	10
21 806 14	25	25	36	70	10
21 806 17	50	32	46	87	10
21 806 24	100	45	56	100	10
21 806 29	150	45	62	105	10
21 806 36	250	45	70	138	10
21 806 44	500	45	86	176	10
21 806 51	750	45	95	203	10
21 806 54	1,000	45	101	225	10
21 806 63	2,000	45	136	260	10
21 806 69	3,500	45	160	295	1
21 806 73	5,000	45	182	330	1
21 806 86	10,000	45	227	410	1
21 806 88	15,000	45	268	445	1
21 806 91	20,000	45	288	505	1

