



single channel PiPettes

## ecopipette Innovation

The **ecopipette** is packed using recyclable, biodegradable materials and constructed from some of the most renewable durable resources available, to ensure an advanced life.

Capp ecotrade™ program: Return any brand of old pipettes in exchange for discounted new ecopipettes™. At Capp ApS, we will ensure an environmentally safe waste discharge of all old pipettes. Ask your favourite authorized Capp distributor for more information about our ecotrade™ program.

### Capp **ecopipette**™ Unique Advantages:

- Lockable volume controller requiring 50% fewer turns
- Display is fully visible during hand operations
- Reduced force requirement for aspiration and blow-out
- The only pipette with fixed and variable exchangeable counter mechanism
- Large pad for ease of tip ejection
- Fully autoclavable pipette body (no disassembly required)
- Double o-ring for optimal tip seal
- Compatible with all major pipette tip brands
- Robust design includes metal tip base
- Lowest incidence of breakage resulting from accidental dropping or tips mounting





# Ordering Info

 **CAPP APS**  
 LANDBRUGSVEJ 10  
 DK-5260 ODENSE  
 DENMARK

TEL. (+45) 6613 6140  
 FAX. (+45) 6613 2770

MAIL: INFO@CAPP.DK  
 WEB: WWW.CAPP.DK

## ecopipette variable volume

Cat No.	Description	Inaccuracy%	Imprecision%	Color
C02-1S	Capp pipette, variable vol. 0.2-2 ul, ecopipette	12.00 / 2.00	6.00 / 1.00	
C02-1	Capp pipette, variable vol. 0.2-2 ul, ecopipette (length as C10-1)	12.00 / 2.00	6.00 / 1.00	
C10-1	Capp pipette, variable vol. 0.5-10 ul, ecopipette	2.50 / 1.00	1.80 / 0.50	
C20-1	Capp pipette, variable vol. 2-20 ul, ecopipette	5.00 / 1.00	1.50 / 0.30	
C50-1	Capp pipette, variable vol. 5-50 ul, ecopipette	2.00 / 1.00	1.20 / 0.30	
C100-1	Capp pipette, variable vol. 10-100 ul, ecopipette	2.50 / 0.70	1.00 / 0.20	
C200-1	Capp pipette, variable vol. 20-200 ul, ecopipette	2.50 / 0.60	0.70 / 0.20	
C1000-1	Capp pipette, variable vol. 100-1000 ul, ecopipette	2.00 / 0.80	0.60 / 0.20	
C5000-1	Capp pipette, variable vol. 1-5 ml, ecopipette	1.20 / 0.60	0.30 / 0.15	

## ecopipette fixed volume (customized volumes available on request)

Cat No.	Description	Inaccuracy%	Imprecision%
C02-1S-001	Capp pipette, Fixed Volume 1 ul, ecopipette	2.50	1.80
C10-1-005	Capp pipette, Fixed Volume 5 ul, ecopipette	1.20	0.80
C10-1-010	Capp pipette, Fixed Volume 10 ul, ecopipette	1.00	0.50
C50-1-020	Capp pipette, Fixed Volume 20 ul, ecopipette	1.00	0.30
C50-1-025	Capp pipette, Fixed Volume 25 ul, ecopipette	0.80	0.35
C100-1-030	Capp pipette, Fixed Volume 30 ul, ecopipette	1.16	0.40
C100-1-050	Capp pipette, Fixed Volume 50 ul, ecopipette	1.00	0.30
C200-1-100	Capp pipette, Fixed Volume 100 ul, ecopipette	0.70	0.20
C200-1-200	Capp pipette, Fixed Volume 200 ul, ecopipette	0.60	0.20
C1000-1-500	Capp pipette, Fixed Volume 500 ul, ecopipette	0.80	0.15
C1000-1-1000	Capp pipette, Fixed Volume 1000 ul, ecopipette	0.80	0.20
C5000-1A-5000	Capp pipette, Fixed Volume 5000 ul, comfort	0.60	0.15
Cxxx-1-yyy	"Other Fixed vol. pipettes (Note: xxx = max pipette volume; yyy = desired volume)"		

## ecopipette Starter Packs

Cat No.	Description
SKE-03-XX	Three pipettes free of choice, stand, racks of tips and timer

**DOMINIQUE DUTSCHER SAS**