






TEL. (+45) 6613 6140 FAX. (+45) 6613 2770

MAIL: INFO@CAPP.DK WEB: WWW.CAPP.DK

Ordering Info



Cat. No.	Description	Inaccuracy %	Imprecision %	Color
B02-1	Capp pipette, variable vol. 0.1-2 ul, Bravo	12.00 / 2.00	6.00 / 1.00	
B10-1	Capp pipette, variable vol. 0.5-10 ul, Bravo	2.50 / 1.00	1.80/ 0.50	HILL COM
B20-1	Capp pipette, variable vol. 2-20 ul, Bravo	5.00 / 1.00	1.50 / 0.30	all index and
B50-1	Capp pipette, variable vol. 5-50 ul, Bravo	2.00 / 1.00	1.20 / 0.30	CHILITED AND
B100-1	Capp pipette, variable vol. 10-100 ul, Bravo	2.50 / 0.70	1.00 / 0.20	Marie Type
B200-1	Capp pipette, variable vol. 20-200 ul, Bravo	2.50 / 0.60	0.70 / 0.20	man y
B1000-1	Capp pipette, variable vol. 100-1000 ul, Bravo	2.00 / 0.80	0.60 / 0.20	1
B5000-1	Capp pipette, variable vol. 1-5 ml, Bravo	1.20 / 0.60	0.30 / 0.15	100 pt
B10000-1	Capp pipette, variable vol. 1-10 ml, Brayo	3.00 / 0.60	0.60 / 0.15	- PA

Capp Bravo fixed volume (customized volumes available on request)

Cat. No.	Description	Inaccuracy %	Imprecision %	
B01-1F	Capp pipette, fixed vol. 1 ul, Bravo	2.50	1.80	
B02-1F	Capp pipette, fixed vol. 2 ul, Bravo	2.00	1.00	
B05-1F	Capp pipette, fixed vol. 5 ul, Bravo	1.20	0.80	
B10-1F	Capp pipette, fixed vol. 10 ul, Bravo	1.00	0.50	
B20-1F	Capp pipette, fixed vol. 20 ul, Bravo	1.00	0.30	
B50-1F	Capp pipette, fixed vol. 50 ul, Bravo	1.00	0.30	
B100-1F	Capp pipette, fixed vol. 100 ul, Bravo	0.70	0.20	
B200-1F	Capp pipette, fixed vol. 200 ul, Bravo	0.60	0.20	
B500-1F	Capp pipette, fixed vol. 500 ul, Bravo	0.80	0.15	
B1000-1F	Capp pipette, fixed vol. 1000 ul, Bravo	0.80	0.20	
B5000-1F	Capp pipette, fixed vol. 5 ml, Bravo	0.60	0.15	
B10000-1F	Capp pipette, fixed vol. 10 ml, Bravo	0.60	0.15	

Capp Bravo Starter Packs

Cat. No. Description

SK-03-BB Three pipettes free of choice, C-05 stand, 3 racks of tips and 1 timer

Pipette stands

Cat. No.	Description
C-05	Pipette stand for up to 5 pipettes

