

The DURAN® GLS 80® wide neck glass bottle is available in three versions – simply choose the most suitable bottle for your requirements.

### THE VERSATILE DURAN® GLS 80® BOTTLE

The wide opening for safe and easy handling of powders and liquids makes the DURAN® GLS 80® bottle a versatile partner for all routine laboratory applications.

Choose from six different bottle sizes, complete with quick-release plastic cap and pouring ring.



Description	Capacity ml	diameter mm	height* mm	Pack unit	Item no.
GLS 80® Laboratory bottle	250	95	105	10	21 860 36 5
GLS 80® Laboratory bottle	500	101	148	10	11 126 27
GLS 80® Laboratory bottle	1.000	101	218	10	11 127 13
GLS 80® Laboratory bottle	2.000	136	248	10	11 127 15
GLS 80® Laboratory bottle	3.500	160	271	1	21 860 69 5
GLS 80® Laboratory bottle	5.000	182	311	1	11 139 49
GLS 80® Laboratory bottle	10.000	227	395	1	11 139 50
GLS 80® Laboratory bottle	20.000	288	483	1	11 139 51
GLS 80® pouring ring, blue, PP		77	7	10	11 601 66

\*without cap