

## THE PLASTIC-COATED DURAN® GLS 80® PROTECT

Ideal for handling particularly valuable or dangerous materials. The plastic coating (PU) on DURAN® GLS 80® Protect provides effective protection during routine laboratory work and retains the contents if the glass becomes damaged.

- scratch-, leak- and shatter-resistant
- high transparency
- temperature range – 30° C to + 140° C
- autoclavable
- suitable for use in microwaves
- stable and durable

Description	Capacity ml	diameter mm	height* mm	Pack unit	Item no.
GLS 80® Laboratory bottle, protect	250	95	105	10	21 865 36 5
GLS 80® Laboratory bottle, protect	500	101	148	10	11 601 52
GLS 80® Laboratory bottle, protect	1,000	101	218	10	11 601 63
GLS 80® Laboratory bottle, protect	2,000	136	248	10	11 601 64
GLS 80® Laboratory bottle, protect	3,500	160	271	1	21 865 69 5
GLS 80® Laboratory bottle, protect	5,000	182	311	1	11 601 65

\*without  
cap

