

Certificate of Conformity



Article no.	Retrace Code/Batch/Lot no.
293102807	N/A

Article description

DURAN® GL 45 Screw Cap, with 2-hose connector (PP) and EDPM gasket

DURAN® Product Certificate

The articles listed under "article number" are hereby confirmed to have been produced using a quality system conforming to ISO 9001 and comply with the DWK Life Sciences quality standards.

Material of construction

We hereby confirm that the blue and grey plastic components of this item are both made from a polypropylene homopolymer. This product includes a sealing ring gasket manufactured from white EDPM rubber 4660.

Declaration of compliance with Food Contact and Additive Regulations

EDPM 4660 rubber

The material used for the manufacture meets the requirements stipulated in the following guidelines/standards:

- FDA 21 CFR 177.2600 e
- 3-A Hygiene-Standards, Class III
- Free of Animal Derived Ingredients (ADI)

Phthalate statement

We certify, that during the manufacturing of the above products, we do not use or intentionally incorporate into it, plasticisers such as Adipates, ESBO, or Phthalates (However the use of DEP, DEHP or DIBP in the catalyst system, used to manufacture the polypropylene polymer resin, may result in traces of these phthalates in the product, typically in concentrations below 1 ppm.)

Statement on BSE/TSE

We confirm that the raw material used in the manufacture of the referenced product is not of animal origin. Accordingly, transmission of animal spongiform encephalopathy agents via the mentioned components is extremely improbable. We therefore state that these materials and the product is to be considered safe to use in food, pharmaceutical and cosmetic contact applications with respect to BSE and TSE transmissions.

Certified by: _____

Mariska Jefferson
Quality Manager

Date: 2023-10-01