

Staubschutzkappe GL45
 Zeichnung: 0817201-GZA-0GL45--..

Gewinde:
 DIN 168-GL 45x4, 3 Gang
 Zeichnungs-Nr. 0303803-GZA-01694--..


Inhaltsmarkierung

176,5 ± 1,5 Höhenlehre

400
 300
 200
 100

DURAN PURE

500 ml ±5%

 DURAN

Made in Germany

Farbe: weiß

0373001-SZ1-B


Retrace Code
 00000000

(28)

Weitergabe sowie Vervielfältigung dieses Dokuments, Verwertung und Mitteilung seines Inhalts sind verboten, soweit nicht ausdrücklich gestattet. Zuwiderhandlungen verpflichten zu Schadenersatz. Alle Rechte für den Fall der Patent-, Gebrauchsmuster- oder Geschmacksmustereintragung vorbehalten.

The reproduction, distribution and utilization of this document as well as the communication of its contents to others without express authorization is prohibited. Offenders will be held liable for the payment of damages. All rights reserved in the event of the grant of a patent, utility model or design.

A	DWK-Logo	21-AUG-2017	Weindorf
Index	Aenderung / changes	Datum / date	Name
zul. Abweichung allowed dev.	Werkstoff / material	Stück quantity	Gewicht weight
Oberfläche surface	Benennung / title	Pos. position	
bearb. drawn	Datum / date	Name	
geprüft checked	01.12.2010	Weindorf	
	31.08.2017	Fries	
	DURAN WHEATON KIMBLE	Form-Nr. mould number	1694/500
DURAN Produktions GmbH & Co. KG www.DWK-LifeSciences.com		Zeichnungsnummer / drawing number	0335520-GZA-01694
		Revision	A
		Blatt/sheet	1/1
		Ers. / repl.f.	Ers. / repl.d:
		Schutzvermerk ISO 16016 beachten. / Refer to protection notice ISO 16016.	INV



4 Reihen Noppen

Ansicht X

X

5 ± 2,5

Ø86 ± 1,5

2,6 +1,4
 -0,6

