

Technical Data Sheet #0064

1. Product Identification

- Article Nr.: **J7800A**
- Description: Superfrost® Plus adhesion slides, with ground edges 90°, clipped corners
- Application: microscope slide for in-vitro-diagnostic use

2. Traceability

- LOT Number: 7 digits
- Position: box and outer carton
- shelflife: 12 months (indicated on box and outer carton)

3. Specification

- Usage: ready to use, single use only
- Raw material: extrawhite glass
- Refractive Index: $n (\lambda = 546,07 \text{ nm}): 1,5171$
- Dimensions (mm):
 - length: 75,36 (+/- 0,43)
 - width: 25,10 (+/- 0,46)
 - thickness: 1,0 (+/- 0,05)

4. Quality

- Norms and directives
 - Directive 98/79/CE on In-vitro-diagnostic medical devices
 - Norm DIN EN 980
 - Norm ISO 15223


Gerhard Menzel GmbH
Postfach 3157 · 38021 Braunschweig, Germany

- Quality assurance by checking random samples on
 - cleanliness
 - dimensions (length, width)
 - absence of glass defects
 - tissue adhesion
 - packaging


5. Storage and usage instructions

- Slides and cover glasses are particularly suitable for the microscopic analysis of blood, urine and tissue samples.
- Slides and cover glasses are ready to use and should not be reused.
- Be aware of the possibility of harm when dealing with the slides and cover glasses, and take the corresponding safety measures, e.g. wearing gloves.
- When using the slides and cover glasses in machines, the additional usage instructions from the manufacturer of the machinery should be followed.
- When using dyes or other chemicals, please note the safety advice from the manufacturer.

- Storage:
 - Keep dry.
 - Avoid large variations in temperature during both storage and usage. Cooling of the product can cause condensation and lead to condensed water forming between the glasses.
 - Please note: Use older products in storage first.
 - Protect the bottom of the carton from dampness.
 - Slide storage periods: see „use before“ date.



Quality Manager



Marketing Manager