## QUICK DEFROSTING EASY & CLEAN

## Tzidefrost+

DEFROSTING VACUUM

## PROMOTE MAINTENANCE OF YOUR FREEZERS -20°C / -80°C / -150°C

The essential tool for any laboratory with a fleet of freezers. Promote a more regular defrosting process, that avoid the turning off the freezer and transfering its contents. In doing so, extend the lifespan of your freezers and reduce your electricity consumption by -30% to -50%

Dimensions: L375mm x H53mm x W390mm - Weight: 9kg Airflow: 40001/m - Sound level: 64dba





**CONTACT@IZIBIO.FR** 13/15 RUE TAIBTOUT 75009 PARIS - FRANCE

> +33 (0) 1 48 01 48 14 www.izibio.fr

## Professional device

Technical specifications specifically optimized for the aspiration of ice and snow at -80°C

Scraper

Scrape and break the frost

25L tank

Vacuum the frost and ice

HEPA filter at the air outlet

Prevent any contamination in case of vacuuming a sample or sample residue

Izidefrost -86°C