


BIOplastics

FULL Skirted 96 x 0.1 ml Well Firm Deck plates

Fit most liquid handling platforms and some (q)PCR cyclers

Game Changer, substitutes Hard Frame Plates




In short

- Extreme Firm Deck, rigid and flat
- Full skirted, wrap resistant
- Robotic handling systems compatible
- Non removeable laser-mark labelled
- Substitutes two component Hard-Frame plates
- Can be closed with any BIOplastics cap format

Formats

- 0.1 ml low profile Full Skirted
- White, frosted natural or clear
- Made from optimized (q)PCR PP, identical to BIOplastics unique extruded recipes
- With or without lasermarked barcode and/or QR coded
- Any side barcoded up to 5 sides
- Customized laser-marking on request



Compatible Platforms and cyclers

- Fits and can be handled by Robotics platforms
- SBS Footprint, SLAS/ANSI compliant
- Stackable, non nesting, flat, robust, wrap resistant
- Fits 96 well (q)PCR cyclers like Bio-Rad, Analytik Jena, Biometra, Qiagen and others
- Cycler compatibility list can be found on www.bioplastics.com

(q)PCR perfect

- Substitutes 2 component plates
- Highly reproducible results
- Low CO₂ Footprint, recyclable PP
- Pyrogen, DNA, RNA, DNase and RNase free



Designed and manufactured in The Netherlands, Patented

Firm Deck range covers: non-, semi- and full skirted plates, semi skirted 8 and 12 tube strips, Tear-Off mats all in Low Profile (0.1 ml) and Regular Profile (0.2 ml) formats.

© 2022 BIOplastics BV, All rights reserved /336622B

Bio-Rad, Qiagen, Analytik Jena and Biometra are registered trademarks of mentioned companies

BIOplastics

Ordering Information

96 x 0.1 ml Well Firm Deck plates **FULL SKIRTED**

	Low profile 0.1 ml	Laser marked	color	Pack
96-well Firm Deck plate Full Skirted clear, natural	B70671FD	no	clear natural	25 plates
96-well Firm Deck plate Full Skirted frosted, natural, barcoded south	L70671FDB	yes	clear natural	25 plates
96-well Firm Deck plate Full Skirted frosted, natural, barcoded north	L70671FDBN	yes	clear natural	25 plates
96-well Firm Deck plate Full Skirted- frosted, natural, barcoded all 4 sites	L70671FDB4	yes	clear natural	25 plates
96-well Firm Deck plate Full Skirted- frosted, white	B70679FD	no	white	25 plates
96-well Firm Deck plate Full Skirted frosted, white, barcoded south	L70679FDB	yes	white	25 plates
96-well Firm Deck2 plate semi skirted frosted, white, barcoded north	L70679FDBN	yes	white	25 plates
96-well Firm Deck2 plate semi skirted frosted, white, barcoded all 4 sites	L70679FDB4	yes	white	25 plates
Related products for closure (q)PCR plates				
Optical wide area 8-cap strip, extra robust	B57801			120 strips
Optical wide area 8-cap strip mat 12 x 8 Tear Off	B57651			50 mats
Optiseal™ Optical Seal for qPCR (adhesive)	157300			100 seals
Optical wide area 12-cap strip mat 8 x 12 Tear Off	B58651			50 mats

Samples available on request

Designed and manufactured in The Netherlands, Patented

Firm Deck range covers: non-, semi- and full skirted plates, semi skirted 8 and 12 tube strips, Tear-Off mats all in Low Profile (0.1 ml) and Regular Profile (0.2 ml) formats.

© 2022 BIOplastics BV, All rights reserved /336622B

Bio-Rad, Qiagen, Analytik Jena and Biometra are registered trademarks of mentioned companies