


# BIOplastics

## FULL Skirted 96 x 0.1 ml Well Firm Deck plates

Fit most liquid handling platforms and some (q)PCR cyclers

Game Changer, substitutes Hard Frame Plates




### In short

- Extreme Firm Deck, rigid and flat
- Full skirted, wrap resistant
- Robotic handling systems compatible
- Non removeable laser-mark labelled
- Substitutes two component Hard-Frame plates
- Can be closed with any BIOplastics cap format

### Formats

- 0.1 ml low profile Full Skirted
- White, frosted natural or clear
- Made from optimized (q)PCR PP, identical to BIOplastics unique extruded recipes
- With or without lasermarked barcode and/or QR coded
- Any side barcoded up to 5 sides
- Customized laser-marking on request



### Compatible Platforms and cyclers

- Fits and can be handled by Robotics platforms
- SBS Footprint, SLAS/ANSI compliant
- Stackable, non nesting, flat, robust, wrap resistant
- Fits 96 well (q)PCR cyclers like Bio-Rad, Analytik Jena, Biometra, Qiagen and others
- Cycler compatibility list can be found on [www.bioplastics.com](http://www.bioplastics.com)

### (q)PCR perfect

- Substitutes 2 component plates
- Highly reproducible results
- Low CO<sub>2</sub> Footprint, recyclable PP
- Pyrogen, DNA, RNA, DNase and RNase free



Designed and manufactured in The Netherlands, Patented

Firm Deck range covers: non-, semi- and full skirted plates, semi skirted 8 and 12 tube strips, Tear-Off mats all in Low Profile (0.1 ml) and Regular Profile (0.2 ml) formats.

© 2022 BIOplastics BV, All rights reserved /336622B

Bio-Rad, Qiagen, Analytik Jena and Biometra are registered trademarks of mentioned companies

# BIOplastics

## Ordering Information

### 96 x 0.1 ml Well Firm Deck plates **FULL SKIRTED**

	Low profile 0.1 ml	Laser marked	color	Pack
96-well Firm Deck plate Full Skirted clear, <b>natural</b>	<a href="#">B70671FD</a>	no	clear natural	25 plates
96-well Firm Deck plate Full Skirted frosted, <b>natural, barcoded south</b>	<a href="#">L70671FDB</a>	yes	clear natural	25 plates
96-well Firm Deck plate Full Skirted frosted, <b>natural, barcoded north</b>	<a href="#">L70671FDBN</a>	yes	clear natural	25 plates
96-well Firm Deck plate Full Skirted- frosted, <b>natural, barcoded all 4 sites</b>	<a href="#">L70671FDB4</a>	yes	clear natural	25 plates
96-well Firm Deck plate Full Skirted- frosted, <b>white</b>	<a href="#">B70679FD</a>	no	white	25 plates
96-well Firm Deck plate Full Skirted frosted, <b>white, barcoded south</b>	<a href="#">L70679FDB</a>	yes	white	25 plates
96-well Firm Deck2 plate semi skirted frosted, <b>white, barcoded north</b>	<a href="#">L70679FDBN</a>	yes	white	25 plates
96-well Firm Deck2 plate semi skirted frosted, <b>white, barcoded all 4 sites</b>	<a href="#">L70679FDB4</a>	yes	white	25 plates
Related products for closure (q)PCR plates				
Optical wide area 8-cap strip, extra robust	<a href="#">B57801</a>			120 strips
Optical wide area 8-cap strip mat 12 x 8 Tear Off	<a href="#">B57651</a>			50 mats
Optiseal™ Optical Seal for qPCR (adhesive)	<a href="#">157300</a>			100 seals
Optical wide area 12-cap strip mat 8 x 12 Tear Off	<a href="#">B58651</a>			50 mats

### Samples available on request

Designed and manufactured in The Netherlands, Patented

Firm Deck range covers: non-, semi- and full skirted plates, semi skirted 8 and 12 tube strips, Tear-Off mats all in Low Profile (0.1 ml) and Regular Profile (0.2 ml) formats.

© 2022 BIOplastics BV, All rights reserved /336622B

Bio-Rad, Qiagen, Analytik Jena and Biometra are registered trademarks of mentioned companies