



# BIOplastics

## 96 x 0.1 ml Well semi skirted Firm Deck2 plates

Fits all brands 0.1 ml qPCR cyclers

the *One and Only* which fits all 0.1 ml (LP) cyclers



### In short

- Extreme Firm Deck rigid and flat
- Semi skirted, wrap resistant
- Robotic handling systems compatible
- Non removeable laser-mark labelled
- Substitutes two component Hard-Frame plates
- Can be closed with any BIOplastics cap format

### Formats

- 0.1 ml low profile Semi Skirted
- White, frosted natural or clear
- Made from optimized (q)PCR PP, identical to BIOplastics unique extruded recipes
- With or without laser marked barcode and/or QR coded
- Any side barcoded up to 5 sides
- Customized laser-marking on request

### Compatible (q)PCR cyclers and platforms

- **Fits** almost all 96 well **qPCR cyclers** which accept **0.1 ml** plates (eg **Roche, Thermo, Bio-Rad, Life Techn, Qiagen**)
- Stackable, no nesting, flat, robust, wrap resistant
- Fits all 96 well PCR cyclers with height adjustable lids
- Fits and can be handled by Robotics platforms
- Cycler compatibility list can be found on [www.bioplastics.com](http://www.bioplastics.com)

### (q)PCR perfect

- Substitutes 2 component plates
- Enables “one plate fits all strategy” to limit result variations while reducing validation efforts at the same time
- highly reproducible results
- Low CO<sub>2</sub> Footprint, recyclable PP
- Pyrogen, DNA, RNA, DNase and RNase free



Designed and manufactured in The Netherlands, Patented

Firm Deck range covers: non-, semi- and full skirted plates, semi skirted 8 and 12 tube strips, Tear-Off mats all in Low Profile (0.1 ml) and Regular Profile (0.2 ml) formats.

© 2022 BIOplastics BV, All rights reserved /336622C

Roche, Thermo, Bio-Rad and Qiagen are registered trademarks of mentioned company

# BIOplastics

## Ordering Information

96 x 0.1 ml Well semi skirted Firm Deck2 plates

the *One and Only* which fits all 0.1 ml (LP) cyclers

	Low profile 0.1 ml	Laser marked	color	Pack
96-well Firm Deck2 plate semi skirted frosted, <b>natural</b> * <sup>1,2,3,4</sup>	<a href="#">B4690FD2</a>	no	clear natural	25 plates
96-well Firm Deck2 plate semi skirted frosted, <b>natural, barcoded south</b> * <sup>1,2,3,4</sup>	<a href="#">L4690FD2LB</a>	yes	clear natural	25 plates
96-well Firm Deck2 plate semi skirted frosted, <b>natural, barcoded north</b> * <sup>1,2,3,4</sup>	<a href="#">L4690FD2LBN</a>	yes	clear natural	25 plates
96-well Firm Deck2 plate semi skirted frosted, <b>natural, barcoded all 4 sites</b> * <sup>1,2,3,4</sup>	<a href="#">L4690FD2LB4</a>	yes	clear natural	25 plates
96-well Firm Deck2 plate semi skirted frosted, <b>white</b> * <sup>1,2,3,4</sup>	<a href="#">B4699FD2</a>	no	white	25 plates
96-well Firm Deck2 plate semi skirted frosted, <b>white, barcoded south</b> * <sup>1,2,3,4</sup>	<a href="#">L4699FD2LB</a>	yes	white	25 plates
96-well Firm Deck2 plate semi skirted frosted, <b>white, barcoded north</b> * <sup>1,2,3,4</sup>	<a href="#">L4699FD2LBN</a>	yes	white	25 plates
96-well Firm Deck2 plate semi skirted frosted, <b>white, barcoded all 4 sites</b> * <sup>1,2,3,4</sup>	<a href="#">L4699FD2LB4</a>	yes	white	25 plates
Related products for closure (q)PCR plates				
Optical wide area 8-cap strip, extra robust	<a href="#">B57801</a>			120 strips
Optical wide area 8-cap strip math 12 x 8 Tear Off	<a href="#">B57651</a>			50 mats
Optiseal™ Optical Seal for qPCR , (adhesive)	<a href="#">157300</a>			100 seals
Optical wide area 12-cap strip math 8 x 12 Tear Off	<a href="#">B58651</a>			50 mats

### Samples available on request

\*<sup>1,2,3,4</sup> Fit-range: 1 Roche (480 / 96 NB), 2 Thermo Fast 0.1 ml cyclers, 3 Bio-rad All , 4 Qiagen quant and others

Compatible cycler list can be found on [www.bioplastics.com](http://www.bioplastics.com)

Designed and manufactured in The Netherlands, Patented

Firm Deck range covers: non-, semi- and full skirted plates, semi skirted 8 and 12 tube strips, Tear-Off mats all in Low Profile (0.1 ml) and Regular Profile (0.2 ml) formats.

© 2022 BIOplastics BV, All rights reserved /336622C

Roche, Thermo, Bio-Rad and Qiagen are registered trademarks of mentioned company