

Syringe Pump



SPLab02-Dual channel



SPLab10-10 channels

SPLab02-LU



Model Number

- | | |
|----------------------|----------------------|
| SPLab01(1 channel) | SPLab02(2 channels) |
| SPLab04(4 channels) | SPLab06(6 channels) |
| SPLab08(8 channels) | SPLab10(10 channels) |
| SPLab12(12 channels) | |



Applicable to different syringes


Features

The syringe pump with infusion/withdrawal mode, stainless steel body, and large-screen chromatic LCD is easy to operate. It can install and fix different types of injectors or syringes in the same time, and it is suitable for transferring micro flow rates with high precision (non-medical).

SPLab01-E



32.4 KV High Voltage Test



Typical Application

- | Special for high voltage electrostatic spinning

Features

- | Metal spray box, The upper pressure block is POM material.
- | Antistatic design, Can tolerate 35KV high pressure, Special for electrostatic spinning.
- | Suitable for small flow rate, high accuracy transmission of liquid.

