

## Flow Rates Peristaltic Pump

## Suitable Pump Head

LabN1, LabN3, LabN6

3 years warranty



YZ1515x



YZ2515x



MC Series



DG Series

### Features

- | 3.2" color LCD screen display.
- | Flow rate and motor speed display in same screen.
- | Timing function, time range 0.1s-9999 hours, can be used for simple dispensing function.


### Typical Application

- | Special for university laboratory and research institute.
- | Ion chromatography and titrator
- | Pilot scale and industry production

### Model Number

| LabN1, LabN3, LabN6

### LabN Series External Control Schematic Diagram



Flow rate

Motor speed

Change direction

Start/stop

Speed knob

Full speed

Escape

LabN Series Interface and Keypad

