



ALPHA Box 136 x 136 Cryo and Storage Boxes


Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 136 x 136 mm. Height = inner height.

Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.

| Colour | 1 or 2 | Heights | | | | | | |
|--------|-------------------------|---------|-------|-------|-------|--------|--------|--------|
| | | 25 mm | 32 mm | 50 mm | 75 mm | 100 mm | 110 mm | 130 mm |
| WH | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| BU | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| YE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| GR | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| RE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |

How to determine the item no.:



Example 50 mm height, white, standard:

KA 1 050 WH 1



Grid Dividers for ALPHA Box

| Cells for vials | Heights | | | | | |
|---------------------------|-------------|-------------|-------------|-------------|-------------|-------------|
| | 15 mm | 25 mm | 30 mm | 40 mm | 65 mm | 80 mm |
| | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. |
| 16x16 up to ca. 7.5 mm Ø | RA10151616N | RA10251616N | RA10301616N | RA10401616N | RA10651616N | RA10801616N |
| 15x15 up to ca. 8.0 mm Ø | RA10151515N | RA10251515N | RA10301515N | RA10401515N | RA10651515N | RA10801515N |
| 14x14 up to ca. 8.7 mm Ø | RA10151414N | RA10251414N | RA10301414N | RA10401414N | RA10651414N | RA10801414N |
| 13x13 up to ca. 9.4 mm Ø | RA10151313N | RA10251313N | RA10301313N | RA10401313N | RA10651313N | RA10801313N |
| 12x12 up to ca. 10.2 mm Ø | RA10151212N | RA10251212N | RA10301212N | RA10401212N | RA10651212N | RA10801212N |
| 11x11 up to ca. 11.2 mm Ø | RA10151111N | RA10251111N | RA10301111N | RA10401111N | RA10651111N | RA10801111N |
| 10x10 up to ca. 12.5 mm Ø | RA10151010N | RA10251010N | RA10301010N | RA10401010N | RA10651010N | RA10801010N |
| 9x9 up to ca. 13.9 mm Ø | RA10150909N | RA10250909N | RA10300909N | RA10400909N | RA10650909N | RA10800909N |
| 8x8 up to ca. 15.7 mm Ø | RA10150808N | RA10250808N | RA10300808N | RA10400808N | RA10650808N | RA10800808N |
| 7x7 up to ca. 17.9 mm Ø | RA10150707N | RA10250707N | RA10300707N | RA10400707N | RA10650707N | RA10800707N |
| 6x6 up to ca. 21.0 mm Ø | RA10150606N | RA10250606N | RA10300606N | RA10400606N | RA10650606N | RA10800606N |
| 5x5 up to ca. 25.5 mm Ø | RA10150505N | RA10250505N | RA10300505N | RA10400505N | RA10650505N | RA10800505N |
| 4x4 up to ca. 32.0 mm Ø | RA10150404N | RA10250404N | RA10300404N | RA10400404N | RA10650404N | RA10800404N |
| 3x3 up to ca. 42.9 mm Ø | RA10150303N | RA10250303N | RA10300303N | RA10400303N | RA10650303N | RA10800303N |

When used with item number given in the overview (beginning with "RA1..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: ALPHA Box of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 10x10 cells.

KA1050WH1 + 0301010N = Cat. No. KA1050WH10301010N



ALPHA Box long2 136 x 268 Cryo and Storage Boxes


Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 136 x 268 mm. Height = inner height.

Other heights can be produced on request.
Grid dividers can be ordered separately or can be supplied fitted.


Matching Racks

Please request!

| Colour | 1 or 2 | Heights | | | | | | |
|-----------|-------------------------|---------|-------|-------|-------|--------|--------|--------|
| | | 25 mm | 32 mm | 50 mm | 75 mm | 100 mm | 110 mm | 130 mm |
| white WH | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| blue BU | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| yellow YE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| green GR | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| red RE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |



How to determine the item no.:



Example 50 mm height, white, standard:

KA 2050WH1

Grid Dividers for ALPHA Box long2

| Cells for vials | Heights | | | | | |
|---------------------------|-------------|-------------|-------------|-------------|-------------|-------------|
| | 15 mm | 25 mm | 30 mm | 40 mm | 65 mm | 80 mm |
| | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. |
| 16x32 up to ca. 7.5 mm Ø | RA20151632N | RA20251632N | RA20301632N | RA20401632N | RA20651632N | RA20801632N |
| 15x30 up to ca. 8.0 mm Ø | RA20151530N | RA20251530N | RA20301530N | RA20401530N | RA20651530N | RA20801530N |
| 14x28 up to ca. 8.7 mm Ø | RA20151428N | RA20251428N | RA20301428N | RA20401428N | RA20651428N | RA20801428N |
| 13x26 up to ca. 9.4 mm Ø | RA20151326N | RA20251326N | RA20301326N | RA20401326N | RA20651326N | RA20801326N |
| 12x24 up to ca. 10.2 mm Ø | RA20151224N | RA20251224N | RA20301224N | RA20401224N | RA20651224N | RA20801224N |
| 11x22 up to ca. 11.2 mm Ø | RA20151122N | RA20251122N | RA20301122N | RA20401122N | RA20651122N | RA20801122N |
| 10x20 up to ca. 12.5 mm Ø | RA20151020N | RA20251020N | RA20301020N | RA20401020N | RA20651020N | RA20801020N |
| 9x18 up to ca. 13.9 mm Ø | RA20150918N | RA20250918N | RA20300918N | RA20400918N | RA20650918N | RA20800918N |
| 8x16 up to ca. 15.7 mm Ø | RA20150816N | RA20250816N | RA20300816N | RA20400816N | RA20650816N | RA20800816N |
| 7x14 up to ca. 17.9 mm Ø | RA20150714N | RA20250714N | RA20300714N | RA20400714N | RA20650714N | RA20800714N |
| 6x12 up to ca. 21.0 mm Ø | RA20150612N | RA20250612N | RA20300612N | RA20400612N | RA20650612N | RA20800612N |
| 5x10 up to ca. 25.5 mm Ø | RA20150510N | RA20250510N | RA20300510N | RA20400510N | RA20650510N | RA20800510N |
| 4x8 up to ca. 32.0 mm Ø | RA20150408N | RA20250408N | RA20300408N | RA20400408N | RA20650408N | RA20800408N |
| 3x6 up to ca. 42.9 mm Ø | RA20150306N | RA20250306N | RA20300306N | RA20400306N | RA20650306N | RA20800306N |

When used with item number given in the overview (beginning with "RA2..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: ALPHA Box long2 of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 10x20 cells.

KA2050WH1 + 0301020N = Cat. No. KA2050WH10301020N



ALPHA Box short 136 x 94 Cryo and Storage Boxes

Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 136 x 94 mm. Height = inner height.


Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.


Matching Racks

Please request!

| Colour | 1 or 2 | Heights | | | | | | |
|-----------|-------------------------|---------|-------|-------|-------|--------|--------|--------|
| | | 25 mm | 32 mm | 50 mm | 75 mm | 100 mm | 110 mm | 130 mm |
| white WH | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| blue BU | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| yellow YE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| green GR | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| red RE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |



How to determine the item no.:



Example 50 mm height, white, standard:

KAK050WH1

Grid Dividers for ALPHA Box short

| Cells for vials | Heights | | | | | |
|--------------------------|-------------|-------------|-------------|-------------|-------------|-------------|
| | 15 mm | 25 mm | 30 mm | 40 mm | 65 mm | 80 mm |
| | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. |
| 16x11 up to ca. 7.5 mm Ø | RAK0151611N | RAK0251611N | RAK0301611N | RAK0401611N | RAK0651611N | RAK0801611N |
| 15x10 up to ca. 8.0 mm Ø | RAK0151510N | RAK0251510N | RAK0301510N | RAK0401510N | RAK0651510N | RAK0801510N |
| 14x9 up to ca. 8.7 mm Ø | RAK0151409N | RAK0251409N | RAK0301409N | RAK0401409N | RAK0651409N | RAK0801409N |
| 13x8 up to ca. 9.4 mm Ø | RAK0151308N | RAK0251308N | RAK0301308N | RAK0401308N | RAK0651308N | RAK0801308N |
| 12x8 up to ca. 10.2 mm Ø | RAK0151208N | RAK0251208N | RAK0301208N | RAK0401208N | RAK0651208N | RAK0801208N |
| 11x7 up to ca. 11.2 mm Ø | RAK0151107N | RAK0251107N | RAK0301107N | RAK0401107N | RAK0651107N | RAK0801107N |
| 10x6 up to ca. 12.5 mm Ø | RAK0151006N | RAK0251006N | RAK0301006N | RAK0401006N | RAK0651006N | RAK0801006N |
| 9x6 up to ca. 13.9 mm Ø | RAK0150906N | RAK0250906N | RAK0300906N | RAK0400906N | RAK0650906N | RAK0800906N |
| 8x5 up to ca. 15.7 mm Ø | RAK0150805N | RAK0250805N | RAK0300805N | RAK0400805N | RAK0650805N | RAK0800805N |
| 7x4 up to ca. 17.9 mm Ø | RAK0150704N | RAK0250704N | RAK0300704N | RAK0400704N | RAK0650704N | RAK0800704N |
| 6x3 up to ca. 21.0 mm Ø | RAK0150603N | RAK0250603N | RAK0300603N | RAK0400603N | RAK0650603N | RAK0800603N |
| 5x3 up to ca. 25.5 mm Ø | RAK0150503N | RAK0250503N | RAK0300503N | RAK0400503N | RAK0650503N | RAK0800503N |
| 4x2 up to ca. 32.0 mm Ø | RAK0150402N | RAK0250402N | RAK0300402N | RAK0400402N | RAK0650402N | RAK0800402N |
| 3x2 up to ca. 42.9 mm Ø | RAK0150302N | RAK0250302N | RAK0300302N | RAK0400302N | RAK0650302N | RAK0800302N |

When used with item number given in the overview (beginning with "RAK..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: ALPHA Box short of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 10x6 cells.

KAK050WH1 + 0301006N = Cat. No. KAK050WH10301006N



BRAVO Box 133 x 133 Cryo and Storage Boxes

Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 133 x 133 mm. Height = inner height.


Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.

| Colour | 1 or 2 | Heights | | | | | | |
|--------|-------------------------|---------|-------|--------|-------|--------|--------|--------|
| | | 25 mm | 32 mm | 50 mm* | 75 mm | 100 mm | 110 mm | 130 mm |
| WH | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| BU | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| YE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| GR | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| RE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |

* Produced traditionally to the American design.
(Actual internal height: 45 mm)

How to determine the item no.:



Example 50 mm height, white, standard:

KB1050WH1



Grid Dividers for BRAVO Box


| Cells for vials | Heights | | | | | |
|---------------------------|-------------|-------------|-------------|-------------|-------------|-------------|
| | 15 mm | 25 mm | 30 mm | 40 mm | 65 mm | 80 mm |
| | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. |
| 16x16 up to ca. 7.2 mm Ø | RB10151616N | RB10251616N | RB10301616N | RB10401616N | RB10651616N | RB10801616N |
| 15x15 up to ca. 7.7 mm Ø | RB10151515N | RB10251515N | RB10301515N | RB10401515N | RB10651515N | RB10801515N |
| 14x14 up to ca. 8.3 mm Ø | RB10151414N | RB10251414N | RB10301414N | RB10401414N | RB10651414N | RB10801414N |
| 13x13 up to ca. 9.0 mm Ø | RB10151313N | RB10251313N | RB10301313N | RB10401313N | RB10651313N | RB10801313N |
| 12x12 up to ca. 9.8 mm Ø | RB10151212N | RB10251212N | RB10301212N | RB10401212N | RB10651212N | RB10801212N |
| 11x11 up to ca. 10.8 mm Ø | RB10151111N | RB10251111N | RB10301111N | RB10401111N | RB10651111N | RB10801111N |
| 10x10 up to ca. 12.0 mm Ø | RB10151010N | RB10251010N | RB10301010N | RB10401010N | RB10651010N | RB10801010N |
| 9x9 up to ca. 13.6 mm Ø | RB10150909N | RB10250909N | RB10300909N | RB10400909N | RB10650909N | RB10800909N |
| 8x8 up to ca. 15.0 mm Ø | RB10150808N | RB10250808N | RB10300808N | RB10400808N | RB10650808N | RB10800808N |
| 7x7 up to ca. 17.3 mm Ø | RB10150707N | RB10250707N | RB10300707N | RB10400707N | RB10650707N | RB10800707N |
| 6x6 up to ca. 20.3 mm Ø | RB10150606N | RB10250606N | RB10300606N | RB10400606N | RB10650606N | RB10800606N |
| 5x5 up to ca. 24.5 mm Ø | RB10150505N | RB10250505N | RB10300505N | RB10400505N | RB10650505N | RB10800505N |
| 4x4 up to ca. 30.8 mm Ø | RB10150404N | RB10250404N | RB10300404N | RB10400404N | RB10650404N | RB10800404N |
| 3x3 up to ca. 41.2 mm Ø | RB10150303N | RB10250303N | RB10300303N | RB10400303N | RB10650303N | RB10800303N |

When used with item number given in the overview (beginning with "RB1..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: BRAVO Box of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 9x9 cells.

KB1050WH1 + 0300909N = Cat. No. KB1050WH10300909N



BRAVO Box long2 133 x 257 Cryo and Storage Boxes

Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 133 x 257 mm. Height = inner height.

Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.


Matching Racks

Please request!

| Colour | 1 or 2 | Heights | | | | | | |
|-----------|-------------------------|---------|-------|--------|-------|--------|--------|--------|
| | | 25 mm | 32 mm | 50 mm* | 75 mm | 100 mm | 110 mm | 130 mm |
| white WH | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| blue BU | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| yellow YE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| green GR | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| red RE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |


* Produced traditionally to the American design.
(Actual internal height: 45 mm)

How to determine the item no.:



Example 50 mm height, white, standard:

KB 2 050 WH 1



Grid Dividers for BRAVO Box long2


| Cells for vials | Heights | | | | | |
|---------------------------|-------------|-------------|-------------|-------------|-------------|-------------|
| | 15 mm | 25 mm | 30 mm | 40 mm | 65 mm | 80 mm |
| | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. |
| 16x32 up to ca. 7.2 mm Ø | RB20151632N | RB20251632N | RB20301632N | RB20401632N | RB20651632N | RB20801632N |
| 15x30 up to ca. 7.7 mm Ø | RB20151530N | RB20251530N | RB20301530N | RB20401530N | RB20651530N | RB20801530N |
| 14x28 up to ca. 8.3 mm Ø | RB20151428N | RB20251428N | RB20301428N | RB20401428N | RB20651428N | RB20801428N |
| 13x26 up to ca. 9.0 mm Ø | RB20151326N | RB20251326N | RB20301326N | RB20401326N | RB20651326N | RB20801326N |
| 12x24 up to ca. 9.8 mm Ø | RB20151224N | RB20251224N | RB20301224N | RB20401224N | RB20651224N | RB20801224N |
| 11x22 up to ca. 10.8 mm Ø | RB20151122N | RB20251122N | RB20301122N | RB20401122N | RB20651122N | RB20801122N |
| 10x20 up to ca. 12.0 mm Ø | RB20151020N | RB20251020N | RB20301020N | RB20401020N | RB20651020N | RB20801020N |
| 9x18 up to ca. 13.6 mm Ø | RB20150918N | RB20250918N | RB20300918N | RB20400918N | RB20650918N | RB20800918N |
| 8x16 up to ca. 15.0 mm Ø | RB20150816N | RB20250816N | RB20300816N | RB20400816N | RB20650816N | RB20800816N |
| 7x14 up to ca. 17.3 mm Ø | RB20150714N | RB20250714N | RB20300714N | RB20400714N | RB20650714N | RB20800714N |
| 6x12 up to ca. 20.3 mm Ø | RB20150612N | RB20250612N | RB20300612N | RB20400612N | RB20650612N | RB20800612N |
| 5x10 up to ca. 24.5 mm Ø | RB20150510N | RB20250510N | RB20300510N | RB20400510N | RB20650510N | RB20800510N |
| 4x8 up to ca. 30.8 mm Ø | RB20150408N | RB20250408N | RB20300408N | RB20400408N | RB20650408N | RB20800408N |
| 3x6 up to ca. 41.2 mm Ø | RB20150306N | RB20250306N | RB20300306N | RB20400306N | RB20650306N | RB20800306N |

When used with item number given in the overview (beginning with "RB2..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: BRAVO Box long2 of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 9x18 cells.

KB2050WH1 + 0300918N = Cat. No. KB2050WH10300918N



BRAVO Box short 133 x 92 Cryo and Storage Boxes

Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 133 x 92 mm. Height = inner height.

Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.


Matching Racks

Please request!

| Colour | 1 or 2 | Heights | | | | | | |
|-----------|-------------------------|---------|-------|--------|-------|--------|--------|--------|
| | | 25 mm | 32 mm | 50 mm* | 75 mm | 100 mm | 110 mm | 130 mm |
| white WH | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| blue BU | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| yellow YE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| green GR | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| red RE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |

* Produced traditionally to the American design.
(Actual internal height: 45 mm)

How to determine the item no.:



Example 50 mm height, white, standard:

KBK050WH1

Grid Dividers for BRAVO Box short

| Cells for vials | Heights | | | | | |
|--------------------------|-------------|-------------|-------------|-------------|-------------|-------------|
| | 15 mm | 25 mm | 30 mm | 40 mm | 65 mm | 80 mm |
| | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. |
| 16x12 up to ca. 7.2 mm Ø | RBK0151612N | RBK0251612N | RBK0301612N | RBK0401612N | RBK0651612N | RBK0801612N |
| 15x11 up to ca. 7.7 mm Ø | RBK0151511N | RBK0251511N | RBK0301511N | RBK0401511N | RBK0651511N | RBK0801511N |
| 14x10 up to ca. 8.3 mm Ø | RBK0151410N | RBK0251410N | RBK0301410N | RBK0401410N | RBK0651410N | RBK0801410N |
| 13x9 up to ca. 9.0 mm Ø | RBK0151309N | RBK0251309N | RBK0301309N | RBK0401309N | RBK0651309N | RBK0801309N |
| 12x8 up to ca. 9.8 mm Ø | RBK0151208N | RBK0251208N | RBK0301208N | RBK0401208N | RBK0651208N | RBK0801208N |
| 11x7 up to ca. 10.8 mm Ø | RBK0151107N | RBK0251107N | RBK0301107N | RBK0401107N | RBK0651107N | RBK0801107N |
| 10x6 up to ca. 12.0 mm Ø | RBK0151006N | RBK0251006N | RBK0301006N | RBK0401006N | RBK0651006N | RBK0801006N |
| 9x6 up to ca. 13.6 mm Ø | RBK0150906N | RBK0250906N | RBK0300906N | RBK0400906N | RBK0650906N | RBK0800906N |
| 8x5 up to ca. 15.0 mm Ø | RBK0150805N | RBK0250805N | RBK0300805N | RBK0400805N | RBK0650805N | RBK0800805N |
| 7x4 up to ca. 17.3 mm Ø | RBK0150704N | RBK0250704N | RBK0300704N | RBK0400704N | RBK0650704N | RBK0800704N |
| 6x4 up to ca. 20.3 mm Ø | RBK0150604N | RBK0250604N | RBK0300604N | RBK0400604N | RBK0650604N | RBK0800604N |
| 5x3 up to ca. 24.5 mm Ø | RBK0150503N | RBK0250503N | RBK0300503N | RBK0400503N | RBK0650503N | RBK0800503N |
| 4x2 up to ca. 30.8 mm Ø | RBK0150402N | RBK0250402N | RBK0300402N | RBK0400402N | RBK0650402N | RBK0800402N |
| 3x2 up to ca. 41.2 mm Ø | RBK0150302N | RBK0250302N | RBK0300302N | RBK0400302N | RBK0650302N | RBK0800302N |

When used with item number given in the overview (beginning with "RBK..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: BRAVO Box short of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 9x6 cells.

KBK050WH1 + 0300906N = Cat. No. KBK050WH10300906N



CellBox Cryo and Storage Boxes for 15 ml and 50 ml tubes

Useful for all tubes up to approx. 22 mm (15 ml) or 34 mm (50 ml) diameter and 125 mm maximum length. Standard box water-repellent fiberboard boxes suitable for ultra-low temperatures. Quality >special< with additional water-proof coating made of plastic. With and without inserted grid divider. Dimensions = External dimensions incl. lid.

| | | 15 ml | | | |
|-----------|-------------------------|---|---|--|--|
| | | CellBox Mini 5 x 5 Cells (122x122x128 mm) | CellBox Maxi 6 x 6 Cells (148x148x128 mm) | | |
| Colour | 1 or 2 | Cat. No. | Cat. No. | | |
| WH | 1 standard 2 special | KCMi1WH1-5x5 KCMi1WH2-5x5 | KCMA1WH1-6x6 KCMA1WH2-6x6 | | |
| BU | 1 standard 2 special | KCMi1BU1-5x5 KCMi1BU2-5x5 | KCMA1BU1-6x6 KCMA1BU2-6x6 | | |
| YE | 1 standard 2 special | KCMi1YE1-5x5 KCMi1YE2-5x5 | KCMA1YE1-6x6 KCMA1YE2-6x6 | | |
| GR | 1 standard 2 special | KCMi1GR1-5x5 KCMi1GR2-5x5 | KCMA1GR1-6x6 KCMA1GR2-6x6 | | |
| RE | 1 standard 2 special | KCMi1RE1-5x5 KCMi1RE2-5x5 | KCMA1RE1-6x6 KCMA1RE2-6x6 | | |

| | | 50 ml | | | |
|-----------|-------------------------|---|---|--|--|
| | | CellBox Mini 3 x 3 Cells (122x122x128 mm) | CellBox Maxi 4 x 4 Cells (148x148x128 mm) | | |
| Colour | 1 or 2 | Cat. No. | Cat. No. | | |
| WH | 1 standard 2 special | KCMi1WH1-3x3 KCMi1WH2-3x3 | KCMA1WH1-4x4 KCMA1WH2-4x4 | | |
| BU | 1 standard 2 special | KCMi1BU1-3x3 KCMi1BU2-3x3 | KCMA1BU1-4x4 KCMA1BU2-4x4 | | |
| YE | 1 standard 2 special | KCMi1YE1-3x3 KCMi1YE2-3x3 | KCMA1YE1-4x4 KCMA1YE2-4x4 | | |
| GR | 1 standard 2 special | KCMi1GR1-3x3 KCMi1GR2-3x3 | KCMA1GR1-4x4 KCMA1GR2-4x4 | | |
| RE | 1 standard 2 special | KCMi1RE1-3x3 KCMi1RE2-3x3 | KCMA1RE1-4x4 KCMA1RE2-4x4 | | |

| | | 15 ml | | | |
|-----------|-------------------------|---|---|--|--|
| | | CellBox Mini long 5 x 10 Cells (122x237x128 mm) | CellBox Maxi long 6 x 12 Cells (148x287x128 mm) | | |
| Colour | 1 or 2 | Cat. No. | Cat. No. | | |
| WH | 1 standard 2 special | KCMi2WH1-5x10 KCMi2WH2-5x10 | KCMA2WH1-6x12 KCMA2WH2-6x12 | | |
| BU | 1 standard 2 special | KCMi2BU1-5x10 KCMi2BU2-5x10 | KCMA2BU1-6x12 KCMA2BU2-6x12 | | |
| YE | 1 standard 2 special | KCMi2YE1-5x10 KCMi2YE2-5x10 | KCMA2YE1-6x12 KCMA2YE2-6x12 | | |
| GR | 1 standard 2 special | KCMi2GR1-5x10 KCMi2GR2-5x10 | KCMA2GR1-6x12 KCMA2GR2-6x12 | | |
| RE | 1 standard 2 special | KCMi2RE1-5x10 KCMi2RE2-5x10 | KCMA2RE1-6x12 KCMA2RE2-6x12 | | |


| | | 50 ml | | | |
|-----------|-------------------------|--|--|--|--|
| | | CellBox Mini long 3 x 6 Cells (122x237x128 mm) | CellBox Maxi long 4 x 8 Cells (148x287x128 mm) | | |
| Colour | 1 or 2 | Cat. No. | Cat. No. | | |
| WH | 1 standard 2 special | KCMi2WH1-3x6 KCMi2WH2-3x6 | KCMA2WH1-4x8 KCMA2WH2-4x8 | | |
| BU | 1 standard 2 special | KCMi2BU1-3x6 KCMi2BU2-3x6 | KCMA2BU1-4x8 KCMA2BU2-4x8 | | |
| YE | 1 standard 2 special | KCMi2YE1-3x6 KCMi2YE2-3x6 | KCMA2YE1-4x8 KCMA2YE2-4x8 | | |
| GR | 1 standard 2 special | KCMi2GR1-3x6 KCMi2GR2-3x6 | KCMA2GR1-4x8 KCMA2GR2-4x8 | | |
| RE | 1 standard 2 special | KCMi2RE1-3x6 KCMi2RE2-3x6 | KCMA2RE1-4x8 KCMA2RE2-4x8 | | |

| | | 15 ml | | | |
|-----------|-------------------------|---|---|--|--|
| | | CellBox Mini without Cells (122x122x128 mm) | CellBox Maxi without Cells (148x148x128 mm) | | |
| Colour | 1 or 2 | Cat. No. | Cat. No. | | |
| WH | 1 standard 2 special | KCMi1WH1 KCMi1WH2 | KCMA1WH1 KCMA1WH2 | | |
| BU | 1 standard 2 special | KCMi1BU1 KCMi1BU2 | KCMA1BU1 KCMA1BU2 | | |
| YE | 1 standard 2 special | KCMi1YE1 KCMi1YE2 | KCMA1YE1 KCMA1YE2 | | |
| GR | 1 standard 2 special | KCMi1GR1 KCMi1GR2 | KCMA1GR1 KCMA1GR2 | | |
| RE | 1 standard 2 special | KCMi1RE1 KCMi1RE2 | KCMA1RE1 KCMA1RE2 | | |

| | | 50 ml | | | |
|-----------|-------------------------|--|--|--|--|
| | | CellBox Mini long without Cells (122x237x128 mm) | CellBox Maxi long without Cells (148x287x128 mm) | | |
| Colour | 1 or 2 | Cat. No. | Cat. No. | | |
| WH | 1 standard 2 special | KCMi2WH1 KCMi2WH2 | KCMA2WH1 KCMA2WH2 | | |
| BU | 1 standard 2 special | KCMi2BU1 KCMi2BU2 | KCMA2BU1 KCMA2BU2 | | |
| YE | 1 standard 2 special | KCMi2YE1 KCMi2YE2 | KCMA2YE1 KCMA2YE2 | | |
| GR | 1 standard 2 special | KCMi2GR1 KCMi2GR2 | KCMA2GR1 KCMA2GR2 | | |
| RE | 1 standard 2 special | KCMi2RE1 KCMi2RE2 | KCMA2RE1 KCMA2RE2 | | |

Grid Dividers for CellBox (only for empty boxes or repeat orders)

| | Cells | Cat. No. | | Cells | Cat. No. | | Cells | Cat. No. | | Cells | Cat. No. |
|-----------|------------|------------|-----------|-----------|------------|-----------|------------|------------|-----------|-----------|------------|
| Mini | 15 ml 5x5 | RC0900505N | Mini | 50 ml 3x3 | RC0900303N | Maxi | 15 ml 6x6 | RC0900606N | Maxi | 50 ml 4x4 | RC0900404N |
| Mini long | 15 ml 5x10 | RC0900510N | Mini long | 50 ml 3x6 | RC0900306N | Maxi long | 15 ml 6x12 | RC0900612N | Maxi long | 50 ml 4x8 | RC0900408N |
| Mini long | 10 ml 8x12 | RC0900812N | | | | | | | | | |



DELTA Box 133 x 133 Cryo and Storage Boxes


Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With extended lid. External dimensions 133 x 133 mm. Height = inner height.

Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.

| Colour | 1 or 2 | Heights | | | | | | |
|--------|-------------------------|---------------|---------------|-------|-------|--------|--------|--------|
| | | 25 mm | 32 mm | 50 mm | 75 mm | 100 mm | 110 mm | 130 mm |
| WH | 1 standard 2 special | see BRAVO Box | see BRAVO Box | 050 | 075 | 100 | 110 | 130 |
| BU | 1 standard 2 special | | | 050 | 075 | 100 | 110 | 130 |
| YE | 1 standard 2 special | | | 050 | 075 | 100 | 110 | 130 |
| GR | 1 standard 2 special | | | 050 | 075 | 100 | 110 | 130 |
| RE | 1 standard 2 special | | | 050 | 075 | 100 | 110 | 130 |

How to determine the item no.:



Example 50 mm height, white, standard:

KD1050WH1



Grid Dividers for DELTA Box

| Cells for vials | Heights | | | | | |
|---------------------------|-------------|-------------|-------------|-------------|-------------|-------------|
| | 15 mm | 25 mm | 30 mm | 40 mm | 65 mm | 80 mm |
| | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. |
| 16x16 up to ca. 7.2 mm Ø | RB10151616N | RB10251616N | RB10301616N | RB10401616N | RB10651616N | RB10801616N |
| 15x15 up to ca. 7.7 mm Ø | RB10151515N | RB10251515N | RB10301515N | RB10401515N | RB10651515N | RB10801515N |
| 14x14 up to ca. 8.3 mm Ø | RB10151414N | RB10251414N | RB10301414N | RB10401414N | RB10651414N | RB10801414N |
| 13x13 up to ca. 9.0 mm Ø | RB10151313N | RB10251313N | RB10301313N | RB10401313N | RB10651313N | RB10801313N |
| 12x12 up to ca. 9.8 mm Ø | RB10151212N | RB10251212N | RB10301212N | RB10401212N | RB10651212N | RB10801212N |
| 11x11 up to ca. 10.8 mm Ø | RB10151111N | RB10251111N | RB10301111N | RB10401111N | RB10651111N | RB10801111N |
| 10x10 up to ca. 12.0 mm Ø | RB10151010N | RB10251010N | RB10301010N | RB10401010N | RB10651010N | RB10801010N |
| 9x9 up to ca. 13.6 mm Ø | RB10150909N | RB10250909N | RB10300909N | RB10400909N | RB10650909N | RB10800909N |
| 8x8 up to ca. 15.0 mm Ø | RB10150808N | RB10250808N | RB10300808N | RB10400808N | RB10650808N | RB10800808N |
| 7x7 up to ca. 17.3 mm Ø | RB10150707N | RB10250707N | RB10300707N | RB10400707N | RB10650707N | RB10800707N |
| 6x6 up to ca. 20.3 mm Ø | RB10150606N | RB10250606N | RB10300606N | RB10400606N | RB10650606N | RB10800606N |
| 5x5 up to ca. 24.5 mm Ø | RB10150505N | RB10250505N | RB10300505N | RB10400505N | RB10650505N | RB10800505N |
| 4x4 up to ca. 30.8 mm Ø | RB10150404N | RB10250404N | RB10300404N | RB10400404N | RB10650404N | RB10800404N |
| 3x3 up to ca. 41.2 mm Ø | RB10150303N | RB10250303N | RB10300303N | RB10400303N | RB10650303N | RB10800303N |

When used with item number given in the overview (beginning with "RB1..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: DELTA Box of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 9x9 cells.

KD1050WH1 + 0300909N = Cat. No. KD1050WH10300909N



MIKE Box 142 x 142 Cryo and Storage Boxes

Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With lid. External dimensions 142 x 142 mm. Height = inner height.

Other heights can be produced on request.


Grid dividers can be ordered separately or can be supplied fitted.

Matching Racks

Please request!

| Colour | 1 or 2 | Heights | | | | | | |
|--------|-------------------------|---------|-------|-------|-------|--------|--------|--------|
| | | 25 mm | 32 mm | 50 mm | 75 mm | 100 mm | 110 mm | 130 mm |
| WH | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| BU | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| YE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| GR | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| RE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |

How to determine the item no.:



Example 50 mm height, white, standard:

KM1050WH1



Grid Dividers for MIKE Box

| Cells for vials | Heights | | | | | |
|---------------------------|-------------|-------------|-------------|-------------|-------------|-------------|
| | 15 mm | 25 mm | 30 mm | 40 mm | 65 mm | 80 mm |
| | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. |
| 16x16 up to ca. 7.7 mm Ø | RM10151616N | RM10251616N | RM10301616N | RM10401616N | RM10651616N | RM10801616N |
| 15x15 up to ca. 8.3 mm Ø | RM10151515N | RM10251515N | RM10301515N | RM10401515N | RM10651515N | RM10801515N |
| 14x14 up to ca. 9.0 mm Ø | RM10151414N | RM10251414N | RM10301414N | RM10401414N | RM10651414N | RM10801414N |
| 13x13 up to ca. 9.7 mm Ø | RM10151313N | RM10251313N | RM10301313N | RM10401313N | RM10651313N | RM10801313N |
| 12x12 up to ca. 10.6 mm Ø | RM10151212N | RM10251212N | RM10301212N | RM10401212N | RM10651212N | RM10801212N |
| 11x11 up to ca. 11.6 mm Ø | RM10151111N | RM10251111N | RM10301111N | RM10401111N | RM10651111N | RM10801111N |
| 10x10 up to ca. 12.8 mm Ø | RM10151010N | RM10251010N | RM10301010N | RM10401010N | RM10651010N | RM10801010N |
| 9x9 up to ca. 14.3 mm Ø | RM10150909N | RM10250909N | RM10300909N | RM10400909N | RM10650909N | RM10800909N |
| 8x8 up to ca. 16.2 mm Ø | RM10150808N | RM10250808N | RM10300808N | RM10400808N | RM10650808N | RM10800808N |
| 7x7 up to ca. 18.6 mm Ø | RM10150707N | RM10250707N | RM10300707N | RM10400707N | RM10650707N | RM10800707N |
| 6x6 up to ca. 21.8 mm Ø | RM10150606N | RM10250606N | RM10300606N | RM10400606N | RM10650606N | RM10800606N |
| 5x5 up to ca. 26.3 mm Ø | RM10150505N | RM10250505N | RM10300505N | RM10400505N | RM10650505N | RM10800505N |
| 4x4 up to ca. 33.0 mm Ø | RM10150404N | RM10250404N | RM10300404N | RM10400404N | RM10650404N | RM10800404N |
| 3x3 up to ca. 44.2 mm Ø | RM10150303N | RM10250303N | RM10300303N | RM10400303N | RM10650303N | RM10800303N |

When used with item number given in the overview (beginning with "RM1..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: MIKE Box of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 10x10 cells.

KM1050WH1 + 0301010N = Cat. No. KM1050WH10301010N



NANU Box 76 x 76 Cryo and Storage Boxes

Matching Racks

See below

Standard box design made of water-repellent fiberboard suitable for ultra-low temperatures. Quality >special< with additional plastic water-proof coating. With extended lid. External dimensions 76 x 76 mm. Height = inner height.


Other heights can be produced on request.

Grid dividers can be ordered separately or can be supplied fitted.

| Colour | 1 or 2 | Heights | | | | | | |
|-----------|-------------------------|---------|-------|-------|-------|--------|--------|--------|
| | | 25 mm | 32 mm | 50 mm | 75 mm | 100 mm | 110 mm | 130 mm |
| white WH | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| blue BU | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| yellow YE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| green GR | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |
| red RE | 1 standard 2 special | 025 | 032 | 050 | 075 | 100 | 110 | 130 |



How to determine the item no.:



Example 50 mm height, white, standard:

KN1050WH1

Grid Dividers for NANU Box

| Cells for vials | Heights | | | | | |
|-------------------------|-------------|-------------|-------------|-------------|-------------|-------------|
| | 15 mm | 25 mm | 30 mm | 40 mm | 65 mm | 80 mm |
| | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. | Cat. No. |
| 7x7 up to ca. 8.0 mm Ø | RN10150707N | RN10250707N | RN10300707N | RN10400707N | RN10650707N | RN10800707N |
| 6x6 up to ca. 10.0 mm Ø | RN10150606N | RN10250606N | RN10300606N | RN10400606N | RN10650606N | RN10800606N |
| 5x5 up to ca. 13.5 mm Ø | RN10150505N | RN10250505N | RN10300505N | RN10400505N | RN10650505N | RN10800505N |
| 4x4 up to ca. 15.0 mm Ø | RN10150404N | RN10250404N | RN10300404N | RN10400404N | RN10650404N | RN10800404N |
| 3x3 up to ca. 21.0 mm Ø | RN10150303N | RN10250303N | RN10300303N | RN10400303N | RN10650303N | RN10800303N |

When used with item number given in the overview (beginning with "RN1..."), grid dividers are always supplied separately (standard supply).

Instead of the customary method of supplying the box and grid divider separately, the box can be supplied with the grid divider already inserted (no surcharge!). For this purpose the code in the item number for the grid divider (highlighted in yellow) must be added to the item number of the box.

Example: NANU Box of height 50 mm, white, standard as well as inserted grid divider of height 30 mm with 5x5 cells.

KN1050WH1 + 0300505N = Cat. No. KN1050WH10300505N

NANU 50 Vertical Racks • 82x82 mm (WxD)

NANU 50 Horizontal Racks • Width 83 mm

>> Useful for NANU 50 and all boxes up to external dimensions 76x76x53 mm (WxDxH). <<

| Cat. No. | Bays | Height (mm) | Cat. No. | Bays | Height (mm) |
|-----------|------|-------------|-----------|------|-------------|
| TN50-04AN | 4 | 226 | TN50-10AN | 10 | 556 |
| TN50-05AN | 5 | 281 | TN50-11AN | 11 | 611 |
| TN50-06AN | 6 | 336 | TN50-12AN | 12 | 666 |
| TN50-07AN | 7 | 391 | TN50-13AN | 13 | 721 |
| TN50-08AN | 8 | 446 | TN50-14AN | 14 | 776 |
| TN50-09AN | 9 | 501 | | | |

| Cat. No. | Bays | Dimensions D x H | Cat. No. | Bays | Dimensions D x H |
|-----------|-------|------------------|-----------|-------|------------------|
| SN50-44AN | 4 x 4 | 329x224 | SN50-65AN | 6 x 5 | 490x279 |
| SN50-45AN | 4 x 5 | 329x279 | SN50-66AN | 6 x 6 | 490x334 |
| SN50-46AN | 4 x 6 | 329x334 | SN50-74AN | 7 x 4 | 571x224 |
| SN50-54AN | 5 x 4 | 410x224 | SN50-75AN | 7 x 5 | 571x279 |
| SN50-55AN | 5 x 5 | 410x279 | SN50-76AN | 7 x 6 | 571x334 |
| SN50-56AN | 5 x 6 | 410x334 | SN50-85AN | 8 x 5 | 651x279 |
| SN50-64AN | 6 x 4 | 490x224 | SN50-86AN | 8 x 6 | 651x334 |

i All vertical racks can be supplied with a safety rod if required. This increases the depth by 5 mm to 87 mm.

Additional vertical and horizontal racks as well as drawer racks for NANU box on request.