


Overall dimensions




Packed weight: 3.5kg
Product weight: 3kg

Packed dimensions;
A 295mm
B 310mm
C 300mm

Product dimensions;
A 185mm
B 200mm
C 200mm

Working Area



| | |
|---|------------------------|
| D | 148mm (at opening) |
| D | 131mm (at tray height) |
| E | 134mm (at opening) |
| E | 117mm (at tray height) |
| F | 150mm |
| G | 132mm |
| H | 150mm |
| I | 187mm |

Vessel capacity with tank lid

| | Diameter | Height | Max. no in bath |
|----------------------|----------|--------|-----------------|
| Duran 21-801: 100ml | 56mm | 105mm | 4 |
| Duran 21-801: 150ml | 62mm | 110mm | 3s |
| Duran 21-801: 250ml | 70mm | 143mm | 2 |
| Duran 21-801: 500ml | 86mm | 181mm | 0 |
| Duran 21-801: 750ml | 95mm | 203mm | 0 |
| Duran 21-801: 1000ml | 101mm | 230mm | 0 |
| Duran 21-801: 2000ml | 136mm | 265mm | 0 |

Vessel capacity without tank lid

| | Diameter | Height | Max. no in bath |
|----------------------|----------|--------|-----------------|
| Duran 21-801: 100ml | 56mm | 105mm | 4 |
| Duran 21-801: 150ml | 62mm | 110mm | 3s |
| Duran 21-801: 250ml | 70mm | 143mm | 2 |
| Duran 21-801: 500ml | 86mm | 181mm | 1 |
| Duran 21-801: 750ml | 95mm | 203mm | 1 |
| Duran 21-801: 1000ml | 101mm | 230mm | 1 |
| Duran 21-801: 2000ml | 136mm | 265mm | 0 |

s - staggered to fit maximum number