


## Overall dimensions



**Packed weight:** 11kg  
**Product weight:** 9kg


### Packed dimensions:

A 445mm  
B 700mm  
C 375mm

### Product dimensions:

A 335mm  
B 590mm  
C 275mm

## Working Area



D	298mm (at opening)
D	281mm (at tray height†)
E	501mm (at opening)
E	485mm (at tray height†)
F	150mm
G	132mm
H	150mm
I	234mm

Vessel capacity with tank lid	Diameter	Height	Max. no in bath
Duran 21-801 : 100ml	56mm	105mm	37s
Duran 21-801 : 150ml	62mm	110mm	28
Duran 21-801 : 250ml	70mm	143mm	24s
Duran 21-801 : 500ml	86mm	181mm	10
Duran 21-801 : 750ml	95mm	203mm	5s
Duran 21-801 : 1000ml	101mm	230mm	0
Duran 21-801 : 2000ml	136mm	265mm	0

Vessel capacity without tank lid	Diameter	Height	Max. no in bath
Duran 21-801 : 100ml	56mm	105mm	37s
Duran 21-801 : 150ml	62mm	110mm	28
Duran 21-801 : 250ml	70mm	143mm	24s
Duran 21-801 : 500ml	86mm	181mm	15
Duran 21-801 : 750ml	95mm	203mm	12s
Duran 21-801 : 1000ml	101mm	230mm	12s
Duran 21-801 : 2000ml	136mm	265mm	6

s - staggered to fit maximum number