

TECHNICAL SPECIFICATIONS

Trade name	Sterilisation indicator tape
Use	Adhesive tape to indicate autoclaving sterilization To be used with individual units (i.e. packs, containers) to demonstrate that the unit has been exposed to the sterilization treatment and to distinguish between processed and unprocessed units. NOT INTENDED FOR USE AS VALIDATION SYSTEM OF STERILIZATION PROCESS.
Material (backing)	Impregnated semicreped paper with indicator strips
Type of adhesive	Natural rubber
Type of sterilization process	Steam 15 min 121°C or 3 min 134°C or 18 min 134°C
Change that occurs after sterilization	Colour of the strips from white to dark brown
Minimum conditions for change	Steam 10 min 121°C or 2 min 134°C

DIMENSIONAL SPECIFICATIONS

Total thickness	mm. 0,15
Tensile strength	10,00%

**GENERAL FEATURES**

CE Mark	No
Conformity	ISO 11140/2005
Product classification	A
Storage conditions after use	Temperature between 15 and 30° C - relative humidity between 30 and 70%
Storage time after use	max. 18 months

PACKAGING

Minimum packaging	8 pcs single wrapped
Box	48 pcs single wrapped

Quality Assurance
Massimiliano Capitanio

DOMINIQUE DUTSHER SAS