




# BOP SI REELS 200m



**Table of combinations double packaging range BOP SI**

Item code	Inner packing (sterile barrier system)	Outer packaging (protective packaging)
02BOPS5020	50mm	80mm
02BOPS8020	80mm	120mm
02BOPS1220	120mm	160mm
02BOPS1620	160mm	220mm
02BOPS2220	220mm	250mm
02BOPS2520	250mm	285mm
02BOPS2820	285mm	320mm
02BOPS3220	320mm	380mm
02BOPS3820	380mm	420mm
02BOPS4220	420mm	-

Echelle 1/3

Indice	Modified by	Date	Draw by	Checked by		
0	Creation	13/03/20	V.Lepesteur	LB	BF	
						Designation
						REEL BOP SI

**\*Information given for reference purposes only. Depending on the nature of the DM to be packaged, the protective packaging may require a size larger than what is indicated in the correspondance below.**