


BOP SI BAGS - EXPORT (SINGLE PACKING)



*Information given for reference purposes only.
Depending on the nature of the DM to be packaged,
the protective packaging may require a size larger
than what is indicated in the correspondance below.


Echelle 1/3					
Indice	Modified by	Date	Draw by	Checked by	
0	Creation	13/03/20	V.Lepesteur	LB BF	
					
					Designation
					BAG BOP SI

Table of combinations double packaging range BOP SI		
Item code	Inner packing (sterile barrier system)	Outer packaging (protective packaging)
90BOPS5020	50x200mm	80x300mm
90BOPS8015	80x150mm	120x250mm
90BOPS8025	80x250mm	120x350mm
90BOPS8030	80x300mm	120x400mm
90BOPS8040	80x400mm	120x500mm
90BOPS1220	120x200mm	160x300mm
90BOPS1225	120x250mm	160x350mm
90BOPS1235	120x350mm	160x450mm
90BOPS1240	120x400mm	160x500mm
90BOPS1250	120x500mm	160x600mm
90BOPS1255	120x550mm	160x650mm

Table of combinations double packaging range BOP SI		
Item code	Inner packing (sterile barrier system)	Outer packaging (protective packaging)
90BOPS1620	160x200mm	220x300mm
90BOPS1630	160x300mm	220x400mm
90BOPS1635	160x350mm	220x450mm
90BOPS1645	160x450mm	220x550mm
90BOPS1650	160x500mm	220x600mm
90BOPS1660	160x600mm	220x700mm
90BOPS1665	160x650mm	220x750mm
90BOP22030	220x300mm	250x400mm
90BOP22040	220x400mm	250x500mm
90BOP22045	220x450mm	250x550mm
90BOP22055	220x550mm	250x650mm

Table of combinations double packaging range BOP SI		
Item code	Inner packing (sterile barrier system)	Outer packaging (protective packaging)
90BOP22060	220x600mm	-
90BOP22070	220x700mm	-
90BOP22075	220x750mm	-
90BOPS2540	250x400mm	285x500mm
90BOPS2550	250x500mm	285x600mm
90BOPS2555	250x550mm	-
90BOPS2565	250x650mm	285x750mm
90BOPS2835	285x350mm	320x450mm
90BOPS2845	285x450mm	320x550mm
90BOPS2850	285x500mm	320x600mm
90BOPS2860	285x600mm	-

Table of combinations double packaging range BOP SI		
Item code	Inner packing (sterile barrier system)	Outer packaging (protective packaging)
90BOPS2875	285x750mm	-
90BOPS3245	320x450mm	380x550mm
90BOPS3255	320x550mm	380x650mm
90BOPS3260	320x600mm	-
90BOPS3850	380x500mm	420x600mm
90BOPS3855	380x550mm	420x650mm
90BOPS3865	380x650mm	-
90BOPS4260	420x600mm	-
90BOPS4265	420x650mm	-