Bellco Glass, Inc. Specification Sheet Product #: 1985-28500

Description: CellectorTM Key Only

- CellectorTM key allows the user to secure and detach mesh screens from CellectorTM Tissue Sieve Homogenizer with ease
- All 304 stainless steel construction
- Pan Dimensions: 1-1/4" OD x 2-1/2" W
- Compatible with 85mL and 130mL collectors

The information contained in this specification sheet we believe is correct to the best of our knowledge. The recommendations and suggestions herein are made without guarantee or representation as to results. We recommend that adequate tests be made in your laboratory or plant to determine if this product meets all of your requirements.



