

Rack for Tubes and Cryo Box



All TPP racks and boxes feature:

- Alpha-numerical inscription to identify your samples
- Heat resistant autoclavable material

Tube Rack with the "click" system:

- Fits any 15 ml and 50 ml centrifuge tube
- Quickly assembled and disassembled
- Stored collapsed saves space
- Without inhibitory burrs on the insertion holes
- Available in 2 sizes

Tissue Culture Tube Rack:

- Fits 8 x tissue culture flat tubes 10,
8 x tissue culture tubes 20
- Transparent top
- Protection rim to protect shifting of tissue culture flat tubes 10

Rack for TubeSpin® Bioreactor 600:


- Fits 3 x tubes
- Reusable
- Stackable
- Safe footprint

Cryo Box:

- Fits 81 x cryo tubes
- Storage of samples at low temperatures $-196\text{ }^{\circ}\text{C}$
- Form stable even in lowest temperatures
- Single handling of lid even after several usages
- Available in 2 sizes

Cryo Rack:


- Fits 40 x cryo tubes
- Ready-mounted rubber feet
- Non-slip properties
- Star-shaped insertion system in the grooves allows one hand operation when used with the star foot TPP cryo tubes





- 1 Rack for TubeSpin® Bioreactor 600.
- 2 The click-system allows easy assembling and disassembling of the racks. A disassembled rack saves storage space.

Type	Product-No.	Capacity	Dimensions l x w x h mm	Material	Qty / Bag	Qty / Case
	99013	3 x 87600	366 x 146 x 70	PP	1	6
	99014	81 x cyro tube 4 ml/5 ml	133 x 133 x 95	PP	1	10
	99015	81 x cyro tube 1.2 ml/2 ml	133 x 133 x 45	PP	1	20
	99016	40 x cyro tube with star-feet	200 x 100 x 25	PC	1	12
	99017	18 x centrifuge tube 15 ml/ 10 x centrifuge tube 50 ml	209 x 98 x 62	PP	1	45
	99018	8 x 91243 or 8 x 91106	200 x 88 x 72	PP	1	28
	99019	30 x centrifuge tube 15 ml/ 20 x centrifuge tube 50 ml	209 x 172 x 62	PP	1	30



- 3 A raised lid design allows the storage of larger cryo tubes.
- 4 The black alpha-numeric scale is placed on the base as well as on the outer side.