



PRODUCT CODE: 255298

Carazzi's Hematoxylin solution (CE-IVD) for clinical diagnostics

Solution for Hematoxylin-Eosin staining, in human and gynaecological samples., In Vitro Diagnostic medical device class A in compliance to the Regulation (EU) 2017/746.

TARIC 3822 00 00 00

PHYSICAL DATA: liquid, Clear, Violet, Miscible with water • D 25/4 1,08 •

BIBLIOGRAPHY: Eduardo Suarez Peregrin., "Manual Técnico de Análisis Clínicos", 9 Ed.,Prieto, 388 (1972) ; J.D. Bancroft, A. Stevens, "Theory and Practice of Histological Techniques", Edinburgh, Scotland •

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l-1,080 kg 1kg-0,926 l

OBSERVATIONS: Filter any sediment present in the solution. •

COMPOSITION:

Hematoxylin	0,1 g
Aluminium Potassium Sulphate 12-hydrate	5 g
Sodium Iodate	0,02 g
Glycerol	20 ml
Water (s.q.m)	100 ml

SPECIFICATIONS:

Identity :
Identity passes test

Maximum limit of impurities

Suitability for hematoxylin - eosin stain passes test

DOMINIQUE DUTSCHER SAS