



Orbitron V™ Platform Rotator

FEATURES

- 3D gyrating motion
- 12" square aluminum tray with flat mat (standard)
- 8" square or 12" square double stacking tray (optional)
- Dimpled rubber mat (optional)
- Supplied with flat mat
- UL/CSA
- CE marked
- 2 year warranty

ACCESSORIES

Catalog No.	Description
260288	Optional 8 inch Tray (Includes Platform, 1 each, flat and dimpled mat)
C 2202208	Dimpled rubber mat 12"
C 2202162	Flat rubber mat 12"
C 2202774	Dimpled rubber mat 8"
C 2202775	Flat rubber mat 8"
260312	Stacking tray 12"

	Orbitron V (115 VAC)	Orbitron V (230 VAC)
Catalog No.	281111	281111-2
Variable speed	5-120 rpm orbits/min	5-120 rpm orbits/min
Variable tilt angle	0° - 12°	0° - 12°
Operating Temperature	Can operate at 4°C to 40°C	Can operate at 4°C to 40°C
Timer	Adjustable up to 12 hours	Adjustable up to 12 hours
Product Dimensions W x D x H	14" x 12" x 11"	14" x 12" x 11"
Load Limit	9 kgs. (20 lbs.)	9 kgs. (20 lbs.)
Power	115 VAC; 50/60Hz, 6 watts, 0.6 amps	230VAC; 50/60Hz, 60 Watts, 0.7 amps CE marked
Shipping Weight	25 lbs.	25 lbs.



Boekel Rocker II™

FEATURES

- Ambient temperature range 4°C to 40°C
- Low voltage power is supplied by wall mount transformer
- Accepts a wide variety of flasks and plates and tubes
- Adjust tilt angle and speed without the need for tools
- Thin power cord easily passes through incubator door gaskets
- Optional stacking tray and dimpled mat (for test tubes) also available
- Supplied with flat non-skid mat
- Power pac is UL, CSA and CE approved
- CE marked
- 2 year warranty

ACCESSORIES

Catalog No.	Description
260355	Stacking Tray (Supplied with rubber non-skid mat)
C2202208	Dimpled Mat

	Rocker II (115 VAC)	Rocker II (230 VAC)
Catalog No.	260350	260350-2
Adjustable Speed	3 - 25 rpm	3 - 25 rpm
Tilt Angle	0° to 20°	0° to 20°
Load Limit	2.7 kg (6 lbs)	2.7 kg (6 lbs)
Operating Temperature	4°C to 40°C	4°C to 40°C
Stacking Tray Weight	0.7 kg (1.5 lbs)	0.7 kg (1.5 lbs)
Stacking Tray Clearance	6.3 cm (2.5")	6.3 cm (2.5")
Product Dimensions w/o platform W x D x H	19.1 cm x 17.8 cm x 12.7 cm (7.5" x 7" x 5")	
Product Dimensions w/platform W x D x H	30.5 cm x 30.5 cm x 17.8 cm (12" x 12" x 7")	
Power	12 V dc, 0.5 A via 115/ 230 V AC adapter CE marked	
Shipping Weight	5.9 kg (13 lbs)	5.9 kg (13 lbs)



855 Pennsylvania Boulevard • Feasterville, PA 19053 USA

800-336-6929 • Phone: 215-396-8200 • Fax: 215-396-8264

E-mail: boekel-info@boekelsci.com • Web site: www.boekelsci.com