Biosafety cabinets

Efficiency, safety and value for routine applications

thermo scientific

MSC-Advantage Biological Safety Cabinets

thermo scientific

Sample and user protection Containment, comfort and convenience

The Thermo Scientific[™] MSC-Advantage[™] Class II Biological Safety Cabinet combines safety, energy efficiency and value for today's laboratory. The cabinet is energy efficient and provides both, personal and product protection from biological hazards and contamination. It can be connected to an external exhaust to provide protection for trace amounts of volatile toxic chemicals that are sometimes used in the cabinet. The cabinet meets EN 12469 for Class II biosafety cabinets.



Containment

- Thermo Scientific[™] SmartFlow[™] technology automates airflow balancing and compensation
- Digital Airflow Verification (DAVe) alerts user to compromised product and operator protection
- H14 HEPA filters remove airborne contaminants
- SmartFlow[™] indicator displays performance status



Comfort

- Low noise level enhances user focus
- 10° sloped front for exceptional comfort and ergonomics
- Spacious work area with comfortable armrests
- Increased light intensity for a brighter work environment



Convenience

- SmartClean front window design to improve cleaningand loading of large equipment
- Intuitive control panel displays safety
 and performance data
- Thermo Scientific[™] SmartPort[™] design organises tubing and cables
- Programmable UV light extends bulb life and saves energy

Basic controls and convenient features



Specifications

| Size/internal width | 0.9m | 1.2m | 1.5m | 1.8m | | |
|--|-------------------|---|-------------------|-------------------|--|--|
| Exterior dimensions WxHxD (mm) | 1000 x 1522 x 798 | 1300 x 1522 x 798 | 1600 x 1522 x 798 | 1900 x 1522 x 798 | | |
| Interior dimensions WxHxD (mm) | 900 x 780 x 630 | 1200 x 780 x 630 | 1500 x 780 x 630 | 1800 x 780 x 630 | | |
| Shipping dimensions WxHxD (mm) | 1110 x 1730 x 925 | 1410 x 1730 x 925 | 1710 x 1730 x 925 | 2010 x 1730 x 925 | | |
| Net weight (Kg) | 170 | 200 | 230 | 280 | | |
| Shipping weight (Kg) | 190 | 225 | 265 | 315 | | |
| Working opening height of front window | 200 mm | | | | | |
| Maximum opening height of front window | 535 mm | | | | | |
| Max load capacity per work surface segment | 25 kg | | | | | |
| Certification to EN12469 | TÜV, LNE | TÜV, LNE | ΤÜV | ΤÜV | | |
| Light power (Lx) | >800 | >850 | >1250 | >1300 | | |
| Receptacles | 2 sockets, ri | right rear wall 4 sockets: 2 on each side rear wall | | | | |
| Noise level (dBA) | <55 | <55 | <59 | <59 | | |
| Energy consumption, operating set point (W) | 150 | 200 | 280 | 340 | | |
| Energy consumption, reduced flow mode (W) | 40 | 40 | 70 | 70 | | |
| Normal heat output in operation (BTU/hr) | 512 | 682 | 955 | 1160 | | |
| Normal heat output in reduced flow mode (BTU/hr) | 136 | 136 | 239 | 239 | | |
| Exhaust/inflow air volume (m3/h) | 292 | 389 | 486 | 583 | | |
| Ducting requirements, thimble connection (m3/h) | 379 | 505 | 632 | 758 | | |





Manually adjustable stand enables cabinet to be set at a comfortable working position between 750-950mm

Configurable work surface improves flexibility of the working area



Side wall service taps are easy to install



Utility basket with hanging bar increases flexibility of the working area and helps improve organisation



SmartPort floor stand hanging shelf to keep area beneath cabinet clear to ensure adequate legroom for a comfortable working position



Direct duct exhaust can be added when working with trace amounts of volatile chemicals

Ordering information and accessories

| Description | Cat. no. | | | | |
|---|-------------------|----------|----------|----------|--|
| | 0.9m | 1.2m | 1.5m | 1.8m | |
| MSC-Advantage Class II Biological Safety Cabinet, 230V, 50/60Hz, EN12469 | 51028225 | 51025411 | 51028226 | 51025413 | |
| Fixed floor stand (height of 680mm) | 50116441 | 50109309 | 50116442 | 50109311 | |
| Castor set - fixed stands only | 50109977 | | | | |
| Manual adjustible floor stand (height of 680-880mm) | 50116443 | 50109312 | 50116444 | 50109313 | |
| Electronic floor stand (height of 680-880mm) | 50117449 | 50109314 | 50117450 | 50109325 | |
| Foot rest bars (Not for electric stands) | 50051982 | 50051983 | 50051984 | 50051985 | |
| Removal of 316 SS work surface | 51901025 | 51900915 | 51901027 | 51900916 | |
| 316 SS worksurface (single piece) | 50076398 | 50082970 | 50122210 | 50085281 | |
| 316 SS worksurface (single piece - lowered) | 50073669 | 50073670 | 50073671 | 50078750 | |
| Replacement ergonomic arm rests (316 SS and gel inserts) | 50156438 | | | | |
| Additional exhaust HEPA filter | | 50109987 | | 50109988 | |
| Additional carbon-enhanced filter | | 50109989 | | 50110305 | |
| Factory-installed UV light option | 51900912 | | | | |
| Vacuum tap G3/8 connection (grey) with long neck | 50044678 | | | | |
| Water tap G3/8 connection (black) with long neck | 50044679 | | | | |
| Non-combustible G3/8 connection (green) with long neck | 50045959 | | | | |
| Combustible G3/8 connection (yellow) with long neck | 50046015 | | | | |
| Solenoid valve for combustible G3/8 (yellow) connection | 50109978 | | | | |
| Lead through for hose connections 10-13mm inner diameter hose | 50059017 | | | | |
| Lead through for hose connections 3mm inner diameter hose | 50076408 | | | | |
| Thimble duct (KTD) | 50109984 50109 | | 50109985 | | |
| Direct duct (KDD) | | 50109981 | | 50109982 | |
| Hanging bar with hooks | 50125639 50109979 | | | | |
| Utility basket for hanging bar kit (hanging bar to be ordered separately) | 50078900 | | | | |
| Ergolign cleanroom saddle stool | ERG-CR-2330 | | | | |
| Floor stand hanging shelf | 50128685 | | | | |
| Replacement SmartPort grommets | 50128681 | | | | |
| Rearwall clips for tubing | 50128703 | | | | |
| Switzerland country configuration | 51900300 | | | | |
| UK country configuration | 51900303 | | | | |
| Italy country configuration | 51900306 | | | | |
| France/ Belgium country configuration | 51900448 | | | | |
| Denmark country configuration | 51900481 | | | | |
| Bulgaria, Czech Republic, Slovakia, Poland country configuration | 51900771 | | | | |
| China country configuration | | 5190 | 0900 | | |

Learn more at thermofisher.com/msc-advantage

thermo scientific

For Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customers' specific uses or applications. © 2023 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **EXT2492 0123**