

Thermo Scientific Savant SpeedVac SPD1030 and SPD2030 Integrated Vacuum Concentrator Systems



Thermo Scientific™ Savant™ SpeedVac™ SPD1030 and SPD2030 Vacuum Concentrators are fully integrated systems that combine concentrator, refrigerated cold trap, diaphragm pump and vacuum gauge in a single, compact unit.



Applications

- DNA/RNA
- Proteomics
- Metabolite research
- Drug and forensics testing
- Environmental sampling

Industries

- Life Science
- Biotech (Molecular Biology and Protein)
- Academic Research
- Pharmaceutical
- Industrial
- Agriculture
- Food

Convenience and flexibility in two capacity ranges

The Savant SpeedVac SPD1030 System is designed for medium capacity (i.e., 120 x 1.5 mL tubes and other formats) while the Savant SpeedVac SPD2030 System accommodates higher capacity (i.e., 200 x 1.5 mL tubes and other formats). A variety of rotors are available for both systems to support sample formats from micro-centrifuge tubes, vials and microplates to large-volume flasks. Both systems fit in a standard laboratory fume hood.

Benefits

- **Proven reliability:** Over 50 years of experience stand behind the Savant SpeedVac family of products. Each unit contains an oil-free PTFE diaphragm vacuum pump, making it virtually maintenance-free.
- **Speed and efficient solvent removal:** Rapid evaporation and efficient solvent removal is enhanced with a radiant cover, while vacuum ramping and pressure controls eliminate sample bumping to maintain sample integrity.
- **Ease of use:** Integrated design and intuitive user interface ensure simplicity and control. The intuitive control panel incorporates a bright display and offers three preset programs and three user-defined programs. Oil-free vacuum source eliminates the need to add or change oil. Data integrity and transfer ability to download data from a USB port.
- **Safety assurance:** SpeedVac Systems' standard glass lids are equipped with a soft-close mechanism that locks securely to help ensure safety for both samples and technicians.

Product Specifications	SpeedVac SPD1030 Integrated System	SpeedVac SPD2030 Integrated System
Physical Properties		
Dimensions (WxDxH)	25 x 26 x 16 in / 64 x 66 x 41 cm	31 x 28 x 19 in / 79 x 71 x 49 cm
Required Ventilation Space	3 in / 8 cm on all sides	3 in / 8 cm on all sides
Shipping Dimensions		
Boxed (WxDxH)	33 x 37 x 31 in / 83 x 94 x 79 cm	33 x 37 x 31 in / 83 x 94 x 79 cm
Product Weight	69 kg / 139 lbs	93 kg / 205 lbs
Shipping Weight (Boxed)	90 kg / 198 lbs	118 kg / 260 lbs
Materials of Construction		
Type of Cover	Radiant glass cover	Radiant glass cover
Chamber	PTFE coated aluminum	PTFE coated aluminum
Chamber Seal	BUNA N-70	BUNA N-70
Fittings	Polypropylene	Polypropylene
Vacuum Pump	2-headed pump with ETFE diaphragms	2-headed pump with ETFE diaphragms
Vacuum Control (Ultimate Vacuum Level)	36 l / 60 Hz, 30 l / 50 Hz <10 Torr (13 mbar, 1.3 kPa)	36 l / 60 Hz, 30 l / 50 Hz <10 Torr (13 mbar, 1.3 kPa)
Operational Requirements		
	SPD1030-115: 115 VAC/60 Hz, 12 Amps SPD1030-230: 230 VAC/50 Hz, 6 Amps	SPD2030-220: 220 VAC/60 Hz, 8 Amps SPD2030-230: 230 VAC/50 Hz, 8 Amps
User Interface/Display		
Temperature Setting	35 - 80 °C in 5 °C increments	35 - 80 °C in 5 °C increments
Two Separate Timer Settings for Heat and Run	1 minute to 9 hours and 59 minutes, continuous	1 minute to 9 hours and 59 minutes, continuous
Vacuum Ramp Rate Setting	1 to 5	1 to 5
Vacuum Level Control	30 to 5.1 torr in 0.1 increments	30 to 5.1 torr in 0.1 increments
Number of Preset Programs	3	3
Number of User-Defined Programs	3	3
Lid Lock Indicator	Yes	Yes
Auto and Manual Run Buttons	Yes	Yes
Radiant Lamp On/Off Button	Yes	Yes
Pre-heat Button	Yes	Yes
Compliance		
	UL listed to U.S. and Canadian requirements. CE marked.	



Small Capacity Rotors for SpeedVac SPD1030 Vacuum Concentrator

Tube Type	Description	No. of Tubes	Working Volumes (mL)	Rotor Part No.
Microcentrifuge Tube	1.5 - 2.0 mL	40	1.2 - 1.6	RH40-11
	1.5 - 2.0 mL	64	1.2 - 1.6	RH64-11
	1.5 - 2.0 mL	120	1.2 - 1.6	RH120-11
Glass and Plastic Tubes	0.5 mL (8 x 29 mm)	100	0.3	RH100-8
	0.4 mL (6 x 50 mm)	40	0.3	RH40-6
	0.4 mL (6 x 50 mm)	100	0.3	RH100-6
	12 x 75 mm, 13 x 100 mm	20	4	RH20-12
	1.5 - 2.0 (12 x 75 mm)	40	4	*RH40-12
	12 x 75 mm	72	4	RH72-12
	13 x 100 mm	10	8	RH20-12
	13 x 100 mm	32	8	RH32-13
	17.5 x 102 mm, 15 mL	8	15	RH8-17.5
	17 x 100 mm, 18 x 100 mm	8	10	RH8-18
	17 x 95 mm, 16 x 100 mm			
	18 x 150 mm	4	24	RH4-18-150
	25 x 105 mm, 30 mL	6	30	RH6-25
	16 x 120 mm, 17 x 120 mm, 15 mL	10	12	RH10-15
	50 mL, 28 x 115 mm	6	40	RH6-50
	12 x 32 mm, 12 x 40 mm	60	2	RH60-12-40
20 x 47 mm, 20 x 58 mm	12	2.4	RH12-20	
Glass Vials	1 dram vial, 15 x 45 mm (4 mL)	24	3	RH24-15
	20 x 60 mm, 20 x 58 mm, 20 x 47 mm	12	4	RH12-20
Mini-Scintillation Vial	18 x 52 mm	24	5.6	RH24-18
	28 x 60 mm	12	16	RH12-28
Pear-Shaped Flask	100 mL	4	80	RH4-100
Microplate	Shallow well plates	2	0.3	RH2MP

*Standard with SPD1030P1 unit

Large Capacity Rotors for SpeedVac SPD2030 Vacuum Concentrator

Tube Type	Description	No. of Tubes	Working Volumes (ml)	Rotor Part No.
Microcentrifuge Tube	1.5 - 2.0 mL	200	1.2 - 1.6	*RH200-12
Glass and Plastic Tubes	12 x 75 mm	200	4	*RH200-12
	13 x 100 mm	118	8	*RH200-12
	16 x 100 mm	96	10	RH96-16-100
	16 x 125 mm (for UMA)	88	12	**RH88-16-125U
	17 x 100 mm, 16 x 100 mm	60	10	RH60-17-100
	18 x 125 mm	48	15	RH48-18-125
	18 x 150 mm, 16 x 125 mm	32	24	RH32-18-150
	15 mL (16 x 120 mm, 17 x 120 mm)	52	12	RH52-15
Centrifuge Tubes	15 mL (17 x 120 mm)	30	12	RH60-17-100
	50 mL (28 x 115 mm), (for UMA)	48	40	**RH48-50U
Glass Vials	1 dram vial, 15 x 45 mm (4 mL)	192	3	RH192-15
	28 x 60 mm, 28 x 50 mm	50	16	RH50-28-60
Flasks	100 mL flask	8	80	RH8-200
	100 mL recovery flask	8	80	RH8-100
Bottles	250 mL centrifuge bottles	8	250	RH8-200
Microplate	Shallow well	12	0.3	MPTR12-210
	Deep well	8	42737	MPTR8-210
Rotor with Swings	4-position rotor holds 4 carriers/4 blocks			**FPR-4A

*Standard with SPD2030P1 unit

**Use with PRO-10, PRO-20 or PRO-50 Upgrades (SUMAX400)

Ordering Information					
SpeedVac Kit	Part No.	Description	SpeedVac Kit	Part No.	Description
SPD1030P1-115	SPD1030-115	Vacuum concentrator	SPD1030P1-230	SPD1030-115	Vacuum concentrator
	RH40-11	Rotor		RH40-11	Rotor
	SCC1	CryoCool™ 1L heat transfer liquid		SCC1	CryoCool 1L heat transfer liquid
	GCF400	Glass condensation flask		GCF400	Glass condensation flask
SPD1030A-115	SPD1030-115	Vacuum concentrator	SPD1030A-230	SPD1030-230	Vacuum concentrator
	GCF400	Glass condensation flask		GCF400	Glass condensation flask
	FC400	Flask cover for GCF400		FC400	Flask cover for GCF400
SPD2030P1-220	SPD2030-220	Vacuum concentrator	SPD2030P1-230	SPD2030-220	Vacuum concentrator
	RH200-12	Rotor		RH200-12	Rotor
	SCC1	CryoCool 1L heat transfer liquid		SCC1	CryoCool 1L heat transfer liquid
	GCF400	Glass condensation flask		GCF400	Glass condensation flask
SPD2030A-220	SPD2030-220	Vacuum concentrator	SPD2030A-230	SPD2030-230	Vacuum concentrator
	GCF400	Glass condensation flask		GCF400	Glass condensation flask
	FC400	Flask cover for GCF400		FC400	Flask cover for GCF400

Rotor Upgrade Kits for SPD2030 Only	
PRO-10	Includes FPR-4A rotor, (4) UPC-1 single-shelf carriers, and SUMAX400 upper magnet assembly
PRO-20	Includes FPR-4A rotor, (4) UPC- 2 two-shelf carriers for microplates, and SUMAX400 upper magnet assembly
PRO-50	Includes FPR-4A rotor, (4) UPC-5 5-shelf carriers for shallow-well microplates, and SUMAX400 upper magnet assembly

Additional Accessories	
Vapor Trap (Included with the System)	
GCF400	Glass condensation flask
FC400	Flask cover for GCF400
145-6012-00	Foam insulating ring
*SCC1	CryoCool heat transfer fluid (1 liter)
Oligo Prep	
ANT100	Oligo prep post-trap assembly
ANS121	Ammonia neutralizing solution (1 bottle) for use with ANT100
Chemical Traps	
SCT120	Chemical trap (order cartridges separately)
DTK120R	Chemical trap kit for radioactivity
*DC120A	Disposable cartridge for neutralizing acid
DC120R	Disposable cartridge for trapping volatile radioactivity
Cart	
CC120/DX	Deluxe convenience cart

*Hazardous Materials -These items require special handling when shipped by air.

Australia +61 39757 4300

Austria +43 1 801 40 0

Belgium +32 9 272 54 82

China +800 810 5118 or

+400 650 5118

France +33 2 2803 2180

Germany national toll free 0800 1 536 376

Germany international +49 6184 90 6321

India toll free 1800 22 8374

India +91 22 6716 2200

Italy +39 02 95059 552

Japan +81 3 5826 1616

Netherlands +31 76 579 55 55

New Zealand +64 9 980 6700

Nordic/Baltic/CIS countries

+358 10 329 2200

Russia +7 812 703 42 15

Spain/Portugal +34 93 223 09 18

Switzerland +41 44 454 12 12

UK/Ireland +44 870 609 9203

USA/Canada +1 866 984 3766

Other Asian countries +852 3107 7600

Countries not listed +49 6184 90 6000

Find out more at thermofisher.com/speedvac

ThermoFisher
SCIENTIFIC