

## Bottles, round, narrow neck, polypropylene



- Translucent, rigid polypropylene bottles with polypropylene screw caps
- Heavy duty and very durable in use
- Autoclavable at 121°C
- Suitable for food contact applications in accordance with EU and USA regulations
- Ideal for a wide range of liquid storage, media preparation and sampling applications
- Batch certificate available to download online
- Good general chemical resistance
- Bulk packs contain uncapped bottles and a pack of caps



Product Code	Capacity (ml)	Height (mm)	O.D. (mm)	Neck I.D. (mm)	Pack qty.
BNP0030PN	30	66	33	13	10
BNP0060PN	60	80	40	13	10
BNP0100PN	100	89	50	18	10
BNP0125PN	125	102	50	18	10
BNP0150PN	150	116	50	18	10
BNP0250P	250	129	62	18	10
BNP0500P	500	167	73	21	10
BNP1000P	1000	213	92	27	5
BNP2000P	2000	240	121	27	5
BNP0030BULK	30	66	32	13	1200
BNP0060BULK	60	80	40	13	800
BNP0100BULK	100	89	50	18	420
BNP0125BULK	125	102	50	18	350
BNP0150BULK	150	116	50	18	330
BNP0250BULK	250	129	60	18	200
BNP0500BULK	500	167	75	21	110
BNP1000BULK	1000	213	90	27	60
BNP2000BULK	2000	240	120	27	30



*When autoclaving bottles always ensure the caps are completely loosened or removed to prevent accidental implosion.*