

### ClearLine syringe filters



The FV syringe filters present an excellent chemical compatibility and resist to organic solvents and strong acids. They are used in pre-filtration of viscous solutions or for filtration of separation of cells media before sterilization

**Features :**

- Pre-filtration of loaded solutions
- Great chemical compatibility

**Specifications :**

- Material : FV – Glass microfiber
- Color code : Natural
- Pore size : 1,0 µm
- Diameter : 33 mm
- Water flow (mL/mn@10psi) : -
- Bubble point (psi) : -
- pH : 1-14
- Sterile : No
- Prefilter : -

**Units :**

100 syringe filters