

### ClearLine syringe filters



Showing a good chemical resistance, MCE syringe filters are efficient in the filtration of proteins, enzymes and other aqueous solutions. The applications are various such as clarification, purification or sterilization of biological fluids.

**Features :**

- Hydrophilic membrane
- Biologically inert

**Specifications :**

- Material : MCE – Mixed cellulose esters
- Color code : Green
- Pore size : 0,45 µm
- Diameter : 25 mm
- Water flow (mL/mn@10psi) : 85-116
- Bubble point (psi) : 25-33
- pH : 4-8
- Sterile : Yes, gamma-ray methods
- Prefilter : -

**Units :**

50 syringe filters