

### ClearLine membranes filters

Naturally hydrophilic, nylon membranes are used for aqueous samples, alkaline or organic filtration of HPLC samples for their chemical resistance to alkaline solutions and solvents.

Membrane : NYL - Polyamide

Pore size (for Ø 47 mm) : 0,45 µm

Water flow @ 10 psi : 9-10 ml/min/cm<sup>2</sup>

Bubble point : 30-32 psi

Diameter : 90 mm

Non sterile

Colour : White, plain

**Units :**

100 membranes