

PHARMED® BPT TUBING / SCHLAUCH

TUBING FOR PERISTALTIC PUMPS AND CELL CULTURE

- Outlasts silicone tubing in peristaltic pumps by up to 30 times
- Can be autoclaved repeatedly
- Heat weldable for sterile access in closed system
- Meets USP Class VI and FDA criteria
- Documented biocompatibility to the ISO 10993 standard
- Density : 0,98 - Shore hardness : 64 - Beige coloured.
- Temperature from -50°C to +135°C (Sterilization : steam).

The american part numbers are indentic to the european Pharmed sizes

Lead time Lieferzeit	European Art. Nr	American Art. Nr	Dimensions / abmessungen						Coil length / Rollen länge (m)
			mm			Inch / Zoll			
			ID	OD / AD	W	ID	OD / AD	W	
2	R650E-26BPT	AY242605	0,50	3,70	1,60	0,02	0,145	0,62	7,5
2	R6501-26BPT	AY242606	0,80	4,00	1,60	1/32.	5/32.	1/16.	7,5
2	R6502-16BPT	AY242002	1,60	3,20	0,80	1/16.	1/8.	1/32.	7,5
2	R6502-26BPT	AY242003	1,60	4,80	1,60	1/16.	3/16.	1/16.	7,5
2	R6503-26BPT	AY242005	2,40	5,60	1,60	3/32.	7/32.	1/16.	7,5
2	R6504-16BPT	AY242006	3,20	4,80	0,80	1/8.	3/16.	1/32.	7,5
2	R6504-26BPT	AY242007	3,20	6,40	1,60	1/8.	1/4.	1/16.	7,5
2	R6506-26BPT	AY242012	4,80	8,00	1,60	3/16.	5/16.	1/16.	7,5
2	R6508-26BPT	AY242017	6,40	9,50	1,60	1/4.	3/8.	1/16.	7,5
2	R6508-36BPT		6,35	11,15	2,40	1/4.	7/16	3,32,	7,5
2	R6508-46BPT	AY242019	6,40	12,70	3,20	1/4.	1/2.	1/8.	7,5
2	R6510-26BPT	AY242022	8,00	11,10	1,60	5/16.	7/16.	1/16.	7,5
2	R6512-26BPT	AY242027	9,50	12,70	1,60	3/8.	1/2.	1/16.	7,5
2	R6512-36BPT		9,50	14,30	2,40	3,8.	9/16,	3,32,	7,5
2	R6512-46BPT	AY242029	9,50	15,90	3,20	3/8.	5/8.	1/8.	7,5
2	R6516-46BPT	AY242038	12,70	19,10	3,20	1/2.	3/4.	1/8.	7,5
2	R6520-46BPT	AY242046	15,90	22,20	3,20	5/8.	7/8.	1/8.	7,5
2	R6524-46BPT	AY242053	19,00	25,40	3,20	3/4.	1.	1/8.	7,5